# EasyOne-CS

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Operator's Manual (ENG)

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Manuale operativo (ITA)

Manual del Usuario (SPA)

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# EasyOne-CS

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## 15 Literatur

- [1] American Thoracic Society. Standardization of Spirometry: 1994 Update, Nov. 11, 1994. Am J Resp Crit Care Med 1995; 152:1107-1136.
- [2] Occupational Health and Safety Administration (OSHA), Pulmonary Function Standards for Cotton Dust, 29 CFR: 1910.1043 Appendix D.
- [3] Social Security Administration Disability (SSD) Guidelines, CFR404: Appendix 1 to Subpart P.
- [4] Ferguson GT, Enright PL, Buist AS, et al. Office spirometry for lung health assessment in adults: A consensus statement from the National Lung Health Education Program. Chest 2000; 117:1146-1161.
- [5] ATS Pulmonary Function Laboratory Management and Procedure Manual, American Thoracic Society, New York, NY 10019.
- [6] Enright PL, Hyatt RE. Office Spirometry. Lea & Febiger, Philadelphia, 1987.
- [7] Hyatt, RE, Scanlon PD, Nakamura M. Interpretation of Pulmonary Function Tests – A Practical Guide. Lippincott – Raven, Philadelphia, 1997.
- [8] American Thoracic Society. Lung Function Testing: Selection of Reference Values and Interpretative Strategies, Am Rev Respir Dis 1991; 144:1202-1218.
- [9] Morris JF, Temple W. Short Report: Spirometric „Lung Age“ Estimation for Motivating Smoking Cessation, Preventive Medicine 14. 655-662 (1985).
- [10] Polgar, Promadhat, Pulmonary Function Testing in Children: Techniques and Standards. W.B. Saunders Co., Philadelphia, 1971.
- [11] M.R. Miller, R. Crapo, J. Hankinson, et al. ATS/ERS Task Force: Standardization of Lung Function Testing. Numbers 1 to 5. Eur Respir J 2005; 26: 153-161, 319-338, 511-522, 720-735, 948-968.
- [12] Global Strategy For The Diagnosis, Management, And Prevention Of Chronic Obstructive Pulmonary Disease, Executive Summary, Updated 2003 (GOLD)
- [13] Hardie et.al., "Risk of over-diagnosis of COPD in asymptomatic elderly never-smokers" Eur Respir J 2002; 20: 1117-1122

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### 17.2 EasyOne-CS Modell 2010 Screen-Connector

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	Referenz	Index	Publikation
Nordamerika	NHANES III (Hankinson)	/REF_P01/	John L. Hankinson, John R. Odencrantz, and Kathleen B. Fedan. Spirometric Reference Values from a Sample of the General U.S. Population. Am J Respir Crit Care Med, Vol 159, p 179-187, 1999
	Knudson_83	/REF_P02/	Knudson, Ronald J, Michael Lebowitz, Holberg Catherine J., Benjamin Burrows. Changes in the Normal Maximal Expiratory Flow-Volume Curve with Aging. American Review of Respiratory Disease, Volume 127, p. 725-734, 1983.
	Knudson_76	/REF_P03/	Knudson, Ronald J, Ronald Slatin, Michael Lebowitz, Benjamin Burrows. The maximal Expiratory Flow-Volume Curve. American Review of Respiratory Disease, Volume 113, p. 587-600, 1976
	Crapo	/REF_P04/	Crapo RO, Morris AH, Gardner RM. Reference spirometric values using techniques and equipment that meets ATS recommendations. Am Rev Respir Dis Volume 123, p.659-664, 1981.
	Morris	/REF_P05/ /REF_P19/	Morris, James F., Koski, Arthur, Lavon Johnson. Spirometric Standards for Healthy Non-Smoking Adults. American Review of Respiratory Disease, Volume 10-3, p. 57-67, 1971 Morris, J.F. West J. Med (1976) 125:110-118
	Hsu	/REF_P06/	Hsu KHK, Bartholomew PH, Thompson V, Hsieh GSJ. Ventilatory Functions of Normal children and Young Adults- Mexican- American, White, Black. I. Spirometry. J Pediatr Volume 95, p. 14-23, 1979.
	Dockery (Harvard)	/REF_P07/	Dockery, D.W. et al., Distribution of Forced Vital Capacity and Forced Expiratory Volume in One Second in Children 6 to 11 Years of Age, American Rev. of Respiratory Disease. Volume 128, p. 405-412, 1983.
	Polgar	/REF_P20/	Polgar, Promadhat, Pulmonary Function Testing in Children: Techniques and Standards. W.B. Saunders Co., Philadelphia, 1971
	Cherniak	/REF_P08/	Cherniak, R.M., and Raber M.B. Normal Standards for Ventilatory Function using an Automated Wedge Spirometer. American Review of Respiratory Disease. Volume 106, p.38-46, 1972
Süd-amerika	Pereira	/REF_P09/	Carlos Alberto de Castro Pereira, Sueli da Penha Barreto, João Geraldo Simões, Francisco W.L. Pereira, José Gerson Gerstler, Joge Nakatani. Valores de referência para a espirometria em uma amostra da população brasileira adulta, Jornal de Pneumologia 18(1):10-22, maio de 1992
Europa	ERS (ECCS, EGKS)	/REF_P10/	P.H. Quanjer. Lung Volumes and Forced Ventilatory Flows. Eur Respir J, Vol 6, Suppl 16, p. 5-40, 1993
	Zapletal	/REF_P11/	A. Zapletal, T. Paul, M. Samanek. Die Bedeutung heutiger Methoden der Lungenfunktionsdiagnostik zur Feststellung einer Obstruktion der Atemwege bei Kindern und Jugendlichen. Z. Erkrank. Atm.-Org., Volume 149, 343-371, 1977.
	Austria (Forche)	/REF_P12/	G. Forche, K. Harnoncourt, E. Stadlober. Neue spirometrische Bezugswerte für Kinder, Jugendliche und Erwachsene. Öst. Ärztztg. 43, 15-16, 1988.
	Sapaldia	/REF_P13/	SAPALDIA team, O Brändli, CH. Schindler, N. Künzli, R. Keller, A.P. Perruchoud. Lung function in healthy never smoking adults: reference values and lower limits of normal of a Swiss population. Thorax 1996; 51:277-283
	Spain (Roca)	/REF_P14/	J. Roca et al. spirometric reference values for a Mediterranean population. Bull Eur Physiopathol Respir, 18:101-102, 1982.
Skandinavien	Hedenström	/REF_P15/	H. Hedenström, P. Malmberg, K. Agarwal. Reference values for lung Function tests in females. Bull. Eur. Physiopathol. Respir. 21, p. 551-557, 1985. H. Hedenström, P. Malmberg, H.V. Fridriksson. Reference values for lung function tests in men. Upsala J. Med. Sci., 91:299-310, 1986
	Gulsvik	/REF_P16/	A. Gulsvik. Spirometri (Korrespondanser). Tidsskr Nor Loegeforen nr. 31, 105:2240-2, 1985.
	Berglund	/REF_P17/	E. Berglund, G. Birath, J. Bjure, G. Grimby, I. Kjellmer, L. Sandqvist, B. Söderholm. Spirometric Studies in Normal Subjects. Acta Medica Scandinavica, Vol. 173, fasc. 2, p. 185-206, 1963.
Australien Asien	Hibbert	/REF_P18/	Marienne E. Hibbert, M App SCI, Anna Lannigan, RN, Louis I. Landau, MD, Peter D. Phelan, MD. Lung Function Values From a Longitudinal Study of Healthy Children Adolescents, Pediatric Pulmonology 7:101-109 (1989)
	Gore, Crockett	/REF_P19/	C.J. Gore, A.J. Crockett, D.G. Pederson, M.L. Booth, A. Bauman, N. Owen. Spirometric standards for healthy adult lifetime nonsmokers in Australia. Eur Respir J., 1995, 8, 773-782
	JRS2001	/REF_P20/	日本人のスパイログラムと動脈血液ガス分圧基準値 日本呼吸器学会肺生理専門委員会 2001年4月

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Address: Technoparkstrasse 1  
CH-8005 Zürich, Switzerland

declares under its sole responsibility, that the product

Product designation: Spirometer

Product name: EasyOne-CS

EasyOne has been classified as Class IIa and is in conformity with the essential requirements and provisions of Council Directive 93/42 EEC,

is in conformity with the following standards transposing harmonized standards:

EN 1041:1998  
EN ISO 14971:2000 and A1 : 2003  
EN 50103:1995  
EN 60601-1: 1990 and A1 :1993 and A2 : 1995  
EN 60601-1-1: 2001  
EN 60601-1-2: 2001  
EN 60601-1-4: 1996 and A1: 1999

is subject to the procedure set out in ISO 9001, ISO 13485, EN 46001 and in Annex 2 of Directive 93/42/EEC under the supervision of Notified Body Number 0120, SGS Yarsley International Certification Services Ltd., 217-221 London Road, Camberley, Surrey, GU15 3EY.

Zurich, March 4<sup>th</sup>, 2005

# EasyOne-CS

## Operator's Manual (ENG)

## Contents

<b>1</b>	<b>Introduction.....</b>	<b>3</b>
<b>2</b>	<b>Warning Information .....</b>	<b>3</b>
<b>3</b>	<b>Intended use.....</b>	<b>4</b>
<b>4</b>	<b>Instrument installation .....</b>	<b>4</b>
4.1	Setting up the device.....	4
4.2	Operation of the keys .....	6
4.3	Setting region, language, date, time, altitude (above sea level) and printer type .....	6
4.4	Installing the PC software .....	7
4.4.1	Installation of the EasyWare PC software .....	7
4.4.2	Installing the Screen-Connector (driver installation) .....	9
4.4.3	Installing EasyWare in a network .....	9
4.4.4	Installing the CardioSoft™ PC software .....	10
4.4.5	Configuring CardioSoft™ .....	10
<b>5</b>	<b>Conducting a spirometry test .....</b>	<b>11</b>
5.1	Preparations .....	11
5.2	Conducting a spirometry test in CardioSoft™ .....	11
5.2.1	Measuring the forced vital capacity (FVC).....	12
5.2.2	Conducting a 'Post' test.....	14
5.2.3	Displaying a detailed procedure report .....	14
5.3	Performing a spirometry test directly with EasyOne .....	15
5.3.1	Measuring the forced vital capacity (FVC).....	15
5.3.2	Adding a spirometry maneuver .....	17
5.3.3	Conducting a 'Post' test.....	17
5.3.4	Saving and retrieving measurements directly in EasyOne.....	17
5.3.5	Quick Test.....	17
5.3.6	Editing patient data.....	18
5.4	Conducting a spirometry test with EasyWare and EasyOne .....	18
<b>6</b>	<b>Further explanation of the spirometry tests .....</b>	<b>18</b>
6.1	Test type 'FVC (Expiration)' .....	19
6.2	Test type 'FVL (Inspiration and Expiration)' .....	19
6.3	Test type 'Slow VC' .....	19
6.4	Test type 'MVV'.....	20
6.5	Test type 'OSHA Cotton Dust'.....	20
6.6	Test type 'Disability' .....	20
6.7	'Post' test.....	20
6.8	Quality messages.....	20
6.9	Interpretation of results.....	21
6.10	Definition of the parameters .....	21
<b>7</b>	<b>EasyWare PC software.....</b>	<b>22</b>
7.1	EasyWare's main window .....	22
7.2	PC database.....	23
7.3	Displaying patient data .....	23
7.4	Entry of patient data .....	23
7.5	Deleting records.....	23
7.6	Test result preview .....	24
7.7	Online test .....	25
7.8	Pediatric animation .....	25
7.9	Trend display .....	25
7.10	Searching for records .....	26
7.11	Printing a test.....	26
7.12	Exporting records .....	26
7.13	Importing records.....	27
7.14	Selecting the EasyOne instrument for the database view .....	27
7.15	Device status .....	27
7.16	Device Configuration .....	27
7.17	Program Preferences .....	27
7.18	Language Selection .....	27

<b>8</b>	<b>Device Configuration.....</b>	<b>28</b>
8.1	Operating modes 'Diagnostic', 'Frontline' and 'NLHEP.'.....	28
8.2	Test Settings .....	28
8.3	System Configuration .....	29
8.4	Print Settings .....	30
<b>9</b>	<b>Quality messages and quality grades .....</b>	<b>31</b>
9.1	Quality messages.....	31
9.2	Quality Grades .....	31
9.3	Best Test Selection.....	32
<b>10</b>	<b>System Interpretation .....</b>	<b>33</b>
10.1	NLHEP System Interpretation .....	33
10.2	GOLD/Hardie System Interpretation.....	34
<b>11</b>	<b>Predicted Values .....</b>	<b>35</b>
11.1	Predicted Values for Adults .....	35
11.2	Predicted Values for Children .....	35
<b>12</b>	<b>Hygiene and instrument maintenance.....</b>	<b>36</b>
<b>13</b>	<b>Calibration Check.....</b>	<b>36</b>
<b>14</b>	<b>Troubleshooting Tips .....</b>	<b>37</b>
<b>15</b>	<b>Bibliography.....</b>	<b>39</b>
<b>16</b>	<b>Accessories .....</b>	<b>39</b>
<b>17</b>	<b>Appendix A: Specifications.....</b>	<b>40</b>
17.1	EasyOne-CS Model 2001 Spirometer .....	40
17.2	EasyOne-CS Model 2010 Screen-Connector.....	40
<b>18</b>	<b>Appendix B: Text file export format.....</b>	<b>41</b>
<b>19</b>	<b>Appendix C: Predicted Normals .....</b>	<b>43</b>
<b>20</b>	<b>Appendix D: Electromagnetic Compatibility (EMC).....</b>	<b>46</b>

## Version history

Version	Description	Date
2030434-002 Ver. A	First edition for EasyOne-CS	27 Sept. 2006
2030434-002 Rev.B	Revised pages: 2, 3, 4, 40, 41, 43 (ECO087884), updated pictures	7 November 2007

## Target reader group

This Operator Manual is intended for specialist medical staff with practical knowledge of medical workflow and applications and familiar with the terminology required for conducting these procedures.

The information in this document may be amended at any time and without prior notice.

Copyright© 2007 by ndd Medizintechnik AG, Switzerland. All rights reserved

EasyOne and the spirette are protected by one or more of the following patents: EP 0597060, EP 0653919, US 5419326, US 5503151, US 5645071 and US5647370.

CardioSoft™ and CASE™ are brands of GE Medical Systems Information Technologies GmbH, a company of the General Electric Company, trading under the name of GE Healthcare.

## 1 Introduction

Thank you for choosing the EasyOne Spirometer.

With EasyOne you have chosen a high quality spirometer that minimizes the need for maintenance due to its unique ultrasound technology. EasyOne does not need calibration and remains consistently accurate over years. The spirette™ breathing tube assures perfectly hygienic conditions for every patient at low cost even if the spirometer is frequently used.

The EasyOne-CS Spirometer contains the Screen-Connector and the related EasyWare PC software, in addition to the CardioSoft™ software used for patient data management. The measurement curves can be displayed in real time on the PC monitor with the EasyOne-CS Screen-Connector.

## 2 Warning Information

The following terms are used as follows in this document:

**Caution:** Possibility of injury or serious damage

**Please note:** Important information for avoiding damage to the instrument or facilitating operation of the instrument

Please note the following information on safe operation of the EasyOne-CS spirometer:

	Important. Note the Operator's Manual.
	Application part, Type BF.
	The product described in this Operator's Manual may not be disposed of as regular, non-separated domestic waste but must be disposed of separately. Please contact an authorized representative of the manufacturer to obtain information on disposal of your unit.
	Catalogue number.
	Serial number.
	E marked per the Medical Device Directive 93/42/EEC of the European Union. The notified body is SGS Yarsley International Certification Services Ltd.
	CSA International classification mark.
	Manufacturer's identification.
	Date of manufacture. The number found under this symbol is the date of manufacture in the YYYY-MM format.
<b>Caution:</b>	The instrument is not suitable for use in the presence of explosive or flammable gases.
<b>Caution:</b>	Connect only computers complying with Standard IEC 60950-1 or IEC 60601-1.

<b>Caution:</b>	For AA batteries, do not attempt to charge, connect improperly, or dispose of in fire as there is possibility of leakage or explosion. Follow manufacturer's recommendation for proper disposal.
-----------------	--

<b>Caution:</b>	Calibration and servicing may be carried out only by GE (General Electric) or ndd staff. Do not open the device.
-----------------	--

<b>Caution:</b>	Pulmonary function tests require maximum effort on the part of the patient and may lead to sensations of dizziness or giddiness.
-----------------	--

<b>Please note:</b>	The unit may be calibrated only by GE Medical Systems Information Technologies After-Sales Servicing Department. Do not open the device.
---------------------	--

<b>Please note:</b>	Please use only long-life batteries (alkaline batteries). Rechargeable batteries may not be used. Remove the batteries from the battery compartment if you do not intend to use the device for long periods. Ensure that the battery compartment cover is properly closed since patient safety may depend on this.
---------------------	--

**Manufacturer:**

ndd Medizintechnik AG  
Technoparkstrasse 1, CH-8005 Zurich, Switzerland  
Tel.: +41 44 445 29 70  
Fax: +41 44 445 25 31  
Web: [www.ndd.ch](http://www.ndd.ch)

**Distributor:**

GE Medical Systems Information Technologies, Inc.  
8200 West Tower Avenue, Milwaukee, WI 53223 USA  
Tel.: +1 414 355 5000  
+800 558 5120 (USA only)  
Fax: +1 414 355 3790

### 3 Intended use

The EasyOne-CS is designed for conducting simple spirometric measurements on adults and children over the age of 4 by general practitioners, specialists, in occupational medicine and in hospitals. The EasyOne-CS spirometer is used together with the spirette™ respiratory tube in order to conduct slow and forced spirometric maneuvers and MVV tests.

## 4 Instrument installation

### 4.1 Setting up the device

The EasyOne Spirometer is supplied with Screen-Connector (USB connection to PC), two alkaline AA batteries, four spirettes, one nose clip, ten nose clip pads, one CD-ROM containing the EasyWare PC software and one CD-ROM containing the CardioSoft™ software.

Install the two AA alkaline batteries (included) in the compartment on the rear of the spirometer, taking care to match the polarity marking on the batteries with the markings inside the battery compartment.

<b>Caution:</b>	Do not attempt to charge or burn the AA batteries used in the instrument. Please follow the manufacturer's instructions on battery disposal.
-----------------	--

<b>Please note:</b>	Please use only long-life batteries (alkaline batteries). Rechargeable batteries may not be used. Remove the batteries from the battery compartment if you do not intend to use the device for long periods. Ensure that the battery compartment cover is properly closed since patient safety may depend on this.
---------------------	--

<b>Please note:</b>	A low battery message will alert you when battery power falls below 10%. Data saved in memory is not lost when battery power is low, or when batteries are removed.
---------------------	---

Install the spirette™ as shown. Be sure to orient the spirette™ so that the arrow on the spirette™ lines up with the arrow on the spirometer.

<b>Caution:</b>	Connect only computers complying with Standard IEC 60950-1 or IEC 60601-1.
-----------------	--

## 4.2 Operation of the keys

- (ON/OFF) This switches EasyOne-CS on or off. Press and hold the key (for at least 2 sec.) until you hear an audible signal.
- (ENTER) This confirms data entry or the selection and moves you to the next entry field.
- (◀) Deletes last character, scrolls to the left or up.
- (▶) Scrolls to the right or down
- (0,ESC) Press the key briefly in order to enter (0), keep the key pressed longer (at least 1 sec.) in order to return to the previous field with (ESC) or abort the operation, press the key briefly twice in order to enter a blank (the key function operates only if letters can be entered).
- (2,abc), etc. Press the key briefly in order to enter the digit "2", press the key briefly in order to enter "A" (the key function operates only if letters can be entered), press the key briefly twice in order to enter "B" (the key function operates only if letters can be entered), if you press the same key quickly several times consecutively, you will scroll first to the upper-case letters, then to the number and then to the lower-case letters, umlauts and special characters can be found on key (1).

**Please note:** The escape key (Esc-0) is particularly helpful and important in unit navigation. The escape function requires the key to be pressed and held momentarily. Escape is useful for moving to previous menus, items, or fields, and escaping a spirometry test. Pressing this key rapidly in fields where letters are possible, such as patient name and report header, allows the entry of a blank space or a zero.

## 4.3 Setting region, language, date, time, altitude (above sea level) and printer type

Press the (ON/OFF) key for at least 2 seconds in order to switch on the instrument. The instrument switches off automatically if no key is pressed for 15 minutes.

You will see a world map on the display when you switch the instrument on for the first time:

Choose your region (North America, South America, Europe, Asia, Africa or Australia) using the arrow keys ◀ ▶ and confirm your entry with ENTER. You are then prompted to choose your language and to set the date, time, contrast, altitude above sea level and the approximate relative humidity at the instrument's location. This data is not preset.

Note that the region can be selected only once. You must use the "Device Configuration / Reset Device" menu to choose a new region.

The spirometer is delivered with the pre-set default settings. Please refer to Chapter 8 of this User Manual for the procedure for changing the settings. Adapting the settings to your needs allows you to get the most out of your EasyOne-CS instrument.

After the above-specified settings have been entered for the first time, all settings apart from the region, can be changed at any time using the DEVICE CONFIGURATION menu item in the main menu. Alternatively, you can change the settings using the PC software.

#### 4.4 Installing the PC software

The minimum requirements for installing the PC software are as follows:


Operating system:	Microsoft Windows 2000 (min. SP4), XP Home and Professional (min. SP2)
Processor:	Pentium $\geq$ 1.6 GHz
Screen resolution:	SVGA 1024 x 768
Hard disk memory:	4 GB (depending on number of procedures), 300 MB free
RAM:	256 MB
Internet Explorer:	Microsoft Internet Explorer 6.0 or higher

(Please refer to the corresponding manual for detailed information for CardioSoft™)


##### 4.4.1 Installation of the EasyWare PC software

You will require Administrator privileges to install the software. Proceed as follows:


- **Do not yet** connect the Screen-Connector to the PC.
- You will require Administrator privileges on the PC in order to install the software.
- Quit all programs and insert the EasyWare CD into your CD ROM drive. The installation routine should start automatically within 30 seconds of inserting the CD. If the installation routine does not start automatically, switch to directory 'Setup' on the CD-ROM and run the Setup.exe program. You will then see the following window:


- Click on "Next >".


- Enter the user name under "User Name" and enter a name for your organization under "Organization". In the bottom part of the window, you can still select whether this application is to be installed for all users ("Anyone who uses this computer") or for you only ("Only for me"). Now click on "Next >".

- Click on "Next >" to confirm the displayed installation path. If you wish to change the displayed installation path, you can do this by clicking on "Change..."

- Confirm the summarized installation information by clicking on "Install".

- Click on "Finish".

- Macromedia Flash Player must be installed on your system in order to display the pediatric animation. Confirm installation by clicking on "Yes".
- Your computer may need to be restarted during this installation procedure under certain circumstances.

- Confirm completion of installation by clicking on "OK". Remove the CD from the drive.

- An EasyWare icon will have been installed on your desktop after installation.
- Launch EasyWare. When EasyWare is launched for the first time, you are prompted to choose the language and enter the software serial number. Enter the software serial number that is provided on the CD-ROM sleeve and then quit EasyWare.

#### 4.4.2 Installing the Screen-Connector (driver installation)

The corresponding driver must be installed in order to install the Screen-Connector. After installation of the EasyWare PC software is complete, connect the Screen-Connector USB connector to a free USB port on your PC.

<b>Caution:</b> Do not touch the contacts of the Screen-Connector.
--

Proceed as follows under Windows XP:

- Message 'Found New Hardware' is displayed.
- The 'Found New Hardware Wizard' is displayed. When prompted whether to make a connection to Windows Update answer "No, not this time". Click on 'Continue'.
- Choose "Install the software automatically (Recommended)". Click on 'Continue'.
- You will see the following message after the driver has been found on the system:


- Accept the proposed driver and choose "Continue Anyway".
- Complete driver installation by choosing "Finish" in the last dialog box. You will now also see a message indicating that the software is ready to use.

Under Windows 2000: Connect the Screen-Connector to a USB port on your PC. Hardware detection is started automatically and the driver is installed. No user entries are necessary.

#### 4.4.3 Installing EasyWare in a network

- Install EasyWare on each PC on which you wish to use EasyWare and set the following configuration.
- You will require Administrator privileges to install the software. During installation, choose the option for allowing all users to use the software. After installation, you will need to start the software once with Administrator privileges. Now enter the software serial number and define the location of the database using Properties in the File menu:

- Enter the complete path and name of the database under 'Database path and name'. You can also search for an existing database with the 'Search' button.
- If you do not yet have a database, enter "EasyWare.mdb" for instance as the name of the database.
- After EasyWare has been installed correctly, you can access an EasyWare database from different PCs at the same time. If you have several EasyOne devices, you can also synchronize several EasyOne devices with the same database at the same time.

Please note: EasyWare is not intended for installation on a CASE™ system. This is why procedures recorded in the network with EasyOne can be displayed only as a summary and not as a detailed report on a CASE™ system.

#### 4.4.4 Installing the CardioSoft™ PC software

You will require Administrator privileges to install the program.

1. Switch on your computer and monitor. Quit all programs (including programs running in the background).
2. Insert the CardioSoft™ CD in the CD-ROM drive. If the installation routine is not started automatically, follow instructions 3 to 5.
3. Click on Start -> Run.
4. Enter the following at the command prompt: X:\Disk1\Setup (X = Name of CD-ROM drive, e.g. "E" or "D").
5. Click on OK.
6. Follow the instructions on screen.
7. Confirm the two proposed directories with *Continue*.
8. Enter the serial number (see CD-ROM). CardioSoft™ is now installed on your computer.
9. Restart Windows.

**Please note:** You will require at least Version V6.5 of CardioSoft™ for operation with EasyWare. The CD-ROM contains the User Manual for CardioSoft™.

#### 4.4.5 Configuring CardioSoft™

You must set the following configuration in order to use CardioSoft™ together with the EasyWare software: Choose 'System Configuration' and then the 'Devices' tab. Make the following settings:

- Path to external program: enter the installation path of EasyWare at this point. Use the button to the right of the entry field for this. Normally, it will be C:\Programs\ndd Medizintechnik\EasyWare\EasyWare.exe.
- Enter the text 'Spirometry' under 'Procedure Name for External Program'.

**Please note:** Do not check the 'Use External Program for Spirometry Devices' checkbox.

- Automatic PDF-Export can be activated in CardioSoft as follows:
  - Select dialog "System Configuration" and click "Export..."
  - Enter the name of your PDF printer driver.
  - Select check box "Generate automatic file name"
    - Select the patient and/or test information to be included in the file name that the system will generate automatically.
  - Select check box "Export after each examination: PDF".

**Please note:** The PDF-Export Directory depends on the settings of the corresponding printer driver.

## 5 Conducting a spirometry test

**Caution:** If operating with the Screen-Connector, always connect EasyOne to the Screen-Connector first before connecting the Screen-Connector to the PC. When disconnecting, always disconnect the USB connector from the PC first before disconnecting the Screen-Connector from EasyOne.

### 5.1 Preparations

Prepare for testing by having the patient loosen tight clothing, remove dentures, and relax. The patient may sit or stand. If standing you may want to perform testing in an area free of sharp table or counter edges, or have a chair handy as there is a slight possibility that the patient could faint during the strenuous spirometry maneuver.

Explain that the purpose of the test is to determine how much air a person's lungs can hold and how quickly that air can be expelled with a forceful, maximal effort. Since the spirometry test requires active participation by the patient it is very important to demonstrate the maneuver for the patient. Emphasize the essential elements of the test:

- filling lungs completely
- sealing lips around the spirette™ so that there are no leaks, taking care not to block its opening with teeth or tongue or bite down excessively
- blasting out as hard and fast as possible
- continue blowing out until the lungs are completely empty


If you are new to spirometry, you should practice testing yourself and others prior to testing patients. You will learn to recognize a poor effort by observing the patient and/or interpreting the Quality Messages displayed by the spirometer after each effort. After a poor effort you must explain what went wrong. Develop enthusiastic coaching techniques to use during the maneuver to maximize your chances of getting quality results with a minimum number of efforts.

**Caution:** Pulmonary function tests require maximum effort on the part of the patient and may lead to sensations of dizziness or giddiness.

### 5.2 Conducting a spirometry test in CardioSoft™

This Chapter describes the procedure for conducting a spirometry test with CardioSoft™, the EasyWare spirometry software and EasyOne-CS.

Launch CardioSoft™. In CardioSoft™, choose a new patient. Select 'Spirometry (external program)' as the type of procedure and click on the 'Select' button. The following dialog box is displayed if EasyOne is not yet switched on or connected to the PC via the Screen-Connector.


Follow the instructions and click on 'OK'.

You will now see a window with the patient data that has been imported directly from CardioSoft™:

If the patient information is incomplete (missing date of birth for instance) or if the patient data lies outside of the permitted range<sup>1</sup>, you will see a corresponding error message.

If necessary, you can add the ethnic group, asthma status and smoker status to the patient data. Confirm the window with OK. You will see the following display:

### 5.2.1 Measuring the forced vital capacity (FVC)


- Insert a spirette™ into the instrument. When doing this, please ensure that the arrow on the spirette™ is lined up with the arrow on the instrument.
- Use the arrow keys to choose the test type, in this case FVC, and confirm your selection with (ENTER).
- Press ENTER when the patient is ready. You will now hear the sensor buzz.
- The instrument prompts you to avoid flow in the spirette™ while it is setting the baseline. It is necessary to block off the spirette™ at one end in order to ensure that the baseline is set precisely. An audible signal sounds when the baseline is set. You will see the prompt "Blast out..." on the screen.
- Hand the instrument to the patient and ask the patient to breath in deeply first, then to insert the spirette™ correctly into his or her mouth, to exhale as firmly and quickly as possible and to

<sup>1</sup> Permitted ranges: age 4 to 120, height 60 to 250 cm, weight 0 to 250 kg.

continue exhaling until all the air has been exhaled. The illustration below shows EasyOne in use (but not using the Screen-Connector).


- The flow volume and the volume-time curve are displayed online on screen during the test. Alternatively, a pediatric animation can be used for FVC tests (see Chapter 7.8).
- The complete flow-volume and volume-time curves are displayed at the end of the spirometry maneuver at the left on screen, and the measured values are shown at the right. An assessment of the measurement and information on how to continue are shown in red beneath the graph with the measured values.

- In order to assess the pulmonary function of the patient, it is necessary to obtain acceptable test quality. The test quality depends on co-operation of the patient and this, in turn, depends on the quality of the physician's instructions. Consequently, EasyOne incorporates an automatic quality control function with prompts to facilitate the physician's job of providing the patient with good instructions. After each maneuver, a message on the screen will inform you as to whether the maneuver was acceptable or not. If not, the message will guide you on how to coach the patient to do better (please see Chapter 6.8 for details):
- Only one of the above prompts is shown after a maneuver. As soon as you see message "Session complete...", you need not conduct further maneuvers. If, even after repeated attempts, it is not possible to obtain an adequate number of good maneuvers, you should take a break, de-

pending on how the patient feels or stop measurement. The results are saved even if the measurement is aborted.

- Now conduct the next test by pressing the (ENTER) key on the EasyOne.
- Comments can be entered even during the test by entering a reviewer name in the text box at the bottom left.
- At least three acceptable spirometry maneuvers must have been performed for a complete spirometry measurement. If this is the case, you will see message "Session Complete! Great Job!". A quality grade from A to F is displayed at the end of the test. It provides information on the overall quality of the test. Please refer to 9.2 for further information on the quality grades.
- You can quit or cancel the test at any time by clicking on "Close".
- Clicking on "Print / Close" ends the test, after which the results are printed out directly.
- After you click on "Close", the data is synchronized between EasyOne and PC. Please wait a moment until this has been done. The following is then displayed on screen:

The screenshot shows a window titled "External Test Summary" with a close button in the top right corner. Below the title bar, it says "Spirometrie[Ext. Prog.] : 20.09.2006 - 17:02:31". The main area contains a table of spirometry results, followed by predicted values, system interpretation, and a confirmation status. At the bottom, there are two buttons: "Detailed Report" and "Cancel".

Parameter	Unit	Value	Pred	%Pred
FEV1	L	3.65	3.63	101
FVC	L	4.08	4.57	89
FEV1/FVC		0.89	0.77	116
MEF75	L	9.92	7.73	128
MEF50	L	6.58	4.74	139
MEF25	L	2.35	1.90	124
MEF2575	L/s	5.23	3.78	138
PEF	L/s	10.56	8.79	120
FET	s	3.48	-	-

Predicted Ref: ERS/Zapletal  
 Syst. Interpret. Ref: NLREP  
 Syst. Interpret.: Normal Spirometry  
 Unconfirmed

The following spirometry tests can be conducted with EasyOne-CS: FVC (expiration), FVL (Inspiration and Expiration), Pre/Post Tests, SVC test and MVV. A detailed description of the test types can be found in Chapter 5.4. The procedure for performing the other test types is equivalent to that for conducting an FVC test.

### 5.2.2 Conducting a 'Post' test


In order to conduct a 'Post' test, choose a patient who has already undergone a spirometry test that day in CardioSoft™. Proceed in precisely the same way as described in the preceding chapter. After choosing the patient, you will then see a window prompting you to specify whether you wish to conduct a 'Post' test or a normal ('Pre') test with this patient. If you answer this prompt with 'Yes', the data of the 'Pre' test will already be shown in the online Test window. Now conduct the 'Post' test on the patient.

If several 'Pre' tests have already been conducted on the same day on a patient, you can choose to what test a 'Post' test is to be added if conducting a 'Post' test.

See Chapter 6.5 for explanatory information on 'Post' test and more detailed information.

### 5.2.3 Displaying a detailed procedure report

Please proceed as follows to display a detailed spirometry procedure report in CardioSoft™. Choose the procedure on a patient in CardioSoft™. Click on the 'View' button. A summary of the external procedure is displayed. Click on the 'Detailed Report' button to view a detailed report. You will see the following display:

You now have the following options:

- **Print:** Choose the 'Print' button on the report display to print out the detailed procedure report. EasyOne offers a number of print options, see Chapter 8.4. The selection of the printer is taken directly from CardioSoft™. If this is not possible, the Windows default printer is used. The header lines on the printout are also taken from CardioSoft™. Choose "System Configuration" for this in CardioSoft™ and enter your data in "Institution Name and Address".
- **Add Comments:** Additional fields for adding comments to the record are shown below the curve display. Enter the name of the reviewer and enter the comments in the large field beneath. After you close with OK, the comments are transferred to CardioSoft™. The comments are also shown on the printout.
- **Delete procedure:** Proceed as follows in order to delete a procedure: you will see a 'Delete' button in the display of the test results at the bottom right. Click on this button and, after you confirm a prompt asking whether you really wish to delete, the test is deleted both from the EasyWare database and from the CardioSoft™ database (after a prompt asking you whether you do really wish to delete).
- **Quit:** Clicking on this button closes the detailed test view and returns you to CardioSoft™.

### 5.3 Performing a spirometry test directly with EasyOne

This Chapter describes the procedure for conducting a spirometry test if using only EasyOne, i.e. EasyOne is not connected to a USB port of the PC via the Screen-Connector. The tests conducted in this way cannot be exported to the CardioSoft™ database.

#### 5.3.1 Measuring the forced vital capacity (FVC)

- Insert a spirette™ into the instrument. When doing this, please ensure that the arrow on the spirette™ is lined up with the arrow on the instrument.
- Choose "Perform Test" in the main menu and then NEW. Confirm with ENTER. The instrument will now allow you to enter the patient data.
- Enter the corresponding patient data line by line. Use the keys as described in Chapter 0. Confirm with ENTER each time.
- Enter the patient data. You then move to the TEST menu. Choose the "FVC-forced Ex test" and confirm with (ENTER).

- Once again briefly prepare the patient for the test. When the patient is ready, press ENTER. You will now hear the sensor buzzing.
- The instrument now prompts you to avoid air flow in the spirette™ since it is setting the baseline. It is advisable to block off the spirette™ on one end in order to ensure that the baseline is set precisely even if the room is draughty. An audible signal will sound when the baseline has been set. You will see prompt "Blast out" on the screen.
- Hand the instrument to the patient. Ask the patient to breathe in deeply, insert the spirette™ correctly into his or her mouth. Now ask the patient to exhale as firmly and as quickly as possible, and continue exhaling until all air has been exhaled.

- The flow-volume curve, the measured values for FEV1 and FVC and the variability of these measured values are displayed after the end of the spirometry maneuver. Confirm this display with (ENTER).
- The next display is used to assess the quality of the measurement. In order to be able to assess the lung function of the patient, it is necessary to achieve a result with acceptable test quality. The test quality depends on the cooperation of the patient and this, in turn, depends on the quality of the instructions given by the examiner. Consequently, EasyOne-CS features automatic quality control with prompts to make it easier for the examiner to give good instructions to the patient. After each spirometry maneuver, a message is displayed on screen, providing information on whether the spirometry maneuver was acceptable and, if it was not, what should be done better. You can find details in Chapter 9.1.
- You can now conduct the next test (CONTINUE), view the data of this test (DATA) and cancel the test (CANCEL) by pressing the arrow keys and (ENTER).
- At least three acceptable spirometry maneuvers must be performed for a complete spirometry measurement. Message "Session Complete! Great Job!" is then displayed. A quality grade from A to F is displayed at the end of the test. It provides information on the overall quality of the test. Please refer to Chapter 9.1 for further information on the quality grades.
- If, even after repeated attempts, it is not possible to obtain an adequate number of good maneuvers, you should take a break, depending on how the patient feels or stop measurement. Even after a break, the measurement is stored and can be printed out under "Print results" in the main menu. You also have the option of adding tests subsequently. Please also read Chapter 5.3.2 for more on this.
- You can view the result on screen (DAT) with the arrow keys. In order to print the result, click on the PRINT button and press (ENTER). Connect the instrument to the PC for this purpose. The report is then printed out.

The following spirometry tests can be conducted with EasyOne-CS: FVC (Expiration), FVL (Inspiration and Expiration), Pre/Post Tests, SVC test and MVV. A detailed description of the test types can be found in Chapter 5.4. The procedure for conducting the other test types is equivalent to that for conducting an FVC test.

### 5.3.2 Adding a spirometry maneuver

If you wish to later add spirometry maneuvers to a test session, for example if the patient was exhausted or if you wanted to test other patients in the meantime, please proceed as follows:

- Select “Perform Test” in the main menu
- Choose the field RECALL and press ENTER
- Scroll through the list of tests until you get to the desired test and press ENTER
- Select the field ADD
- Perform the additional spirometry maneuvers

Please note that it is only possible to add a trial to a previous test that was performed on the same day. No spirometry maneuvers can be added to tests that were generated with CardioSoft™.

### 5.3.3 Conducting a 'Post' test

In order to add a Post test after an FVC or FVL test, click on the POST button on the Result screen.

Proceed as follows if you wish to conduct a Post test from the main menu:

- Select “Perform Test” in the main menu
- Choose the field RECALL and press ENTER
- Scroll through the list of tests until you get to the desired pre-med test and press ENTER
- Select the field POST
- Now proceed as described in Chapter 5.3.1.

### 5.3.4 Saving and retrieving measurements directly in EasyOne

EasyOne-CS saves all test results automatically to the spirometer’s internal memory. No data can be lost even if the batteries are removed from the instrument. If the internal memory of the EasyOne-CS is full (after approx. 700 normal spirometry measurements), the oldest test in the internal memory of the instrument is overwritten<sup>2</sup>.

You can recall saved measurements for the purpose of conducting a new test with the same patient, adding maneuvers, adding a post-test, viewing results, or printing results.

You can add a maneuver or a post-test only on the same day that the original test was performed. Please read Chapter 6.5 und 5.3.2.

In order to add a test to an old measurement, choose “Perform Test” in the main menu and then choose RECALL. Follow the rest of the instructions.

In order to view an old test, choose “View Results” in the main menu and choose the desired test.

### 5.3.5 Quick Test

You have the option to perform a quick test without entering patient data. Select “Perform Test” in the main menu and then QUICK. Choose the desired test using the arrow keys (←) or (→) and perform the test.

<b>Please note:</b> When Quick Test is selected, no comparison to predicted normals are displayed or printed. Predicted normals are only available when age, height and gender are entered.
---

It is possible to enter patient data after having performed a quick test. Proceed as described in Chapter 5.3.6. Once patient data is entered predicted normals will be displayed and printed.

---

<sup>2</sup> No data is deleted or overwritten when synchronizing with the PC. Overwriting tests relates only to test data in the EasyOne-CS’s internal memory.

### 5.3.6 Editing patient data

Patient data exported from Cardiosoft™ to EasyOne cannot be edited. Patient data entered in EasyOne independently of Cardiosoft™ can be edited or supplemented subsequently. In order to do this, go to EDIT DATABASE in the main menu and press (ENTER). Choose the desired measurement with the arrow keys and make the changes.

**Please note:** Editing patient data may influence predicted normal calculation and the system interpretation of the test result. Consequently, recheck the measurement result if height, age, ethnic group or sex are changed.

## 5.4 Conducting a spirometry test with EasyWare and EasyOne

This Chapter describes the procedure for conducting a spirometry test using EasyWare and EasyOne. For this purpose, EasyOne must be operated in so-called Online mode. EasyOne must be connected to a USB port of the PC via the Screen-Connector. The tests conducted in this way cannot be exported to the Cardiosoft™ database.

Please proceed as follows in order to conduct a test:

- Create a new patient record by choosing New Patient in the Edit menu. This step is optional. A new patient record can also be entered subsequently with the EasyOne keyboard (see below).
- Use the menu or the toolbar to switch to Online mode. A window similar to the Test Preview window is displayed on the PC screen.
- The EasyOne can now be operated normally with the instrument's keyboard. Choose "Perform Test / Search" to retrieve the previously entered patient data or choose "Perform Test / New" to enter new patient data. Now conduct a spirometry test as normal.
- When you conduct the test on the patient, flow/volume and volume/time diagrams are shown in real time. The most recent curve is always shown in dark red, and the messages for the current test are shown beneath the numeric data. Alternatively, a pediatric animation can be used for FVC tests (see Chapter 7.8).

Remarks:

- Online mode can be quit at any time with the OK button.
- Online mode can be used with all test types and also with Post test.
- EasyOne synchronizes again automatically with EasyWare when you quit the Online mode.

## 6 Further explanation of the spirometry tests

When you enter the patient data or select an existing patient, you will see the Test menu with the following selection options:

- FVC (expiration)
- FVL (inspiration and expiration)
- MVV
- Slow VC

It is also possible to add a "Post" test to a FVC or FVL Test.

The various measurement methods are outlined below. Good co-operation on the part of the patient is essential with all methods. Consequently, you should explain to the patient clearly beforehand what he or she has to do and motivate the patient to co-operate. Choose the required measurement method with keys (>) or (<) and confirm with ENTER.

## 6.1 Test type 'FVC (Expiration)'

This is the most commonly used spirometric measurement. Prepare the patient as described in Chapter 5.1. before you start the test. Then proceed as follows.

- Insert a spirette™ into the instrument. When doing this, please ensure that the arrow on the spirette™ is lined up with the arrow on the instrument.
- Press ENTER when the patient is ready. You will now hear the sensor buzz.
- The instrument prompts you to avoid flow in the spirette™ while it is setting the baseline. It is necessary to block off the spirette™ at one end in order to ensure that the baseline is set precisely. An audible signal sounds when the baseline is set. You will see the prompt "Blast out..." on the screen.
- Hand the instrument to the patient and ask the patient to breath in deeply first, then to insert the spirette™ correctly into his or her mouth, to exhale as firmly and quickly as possible and to continue exhaling until all the air has been exhaled.
- A message indicating whether the spirometry maneuver was acceptable is displayed at the end of the spirometry maneuver. At least three acceptable, reproducible spirometry maneuvers must be performed until message "Session Complete! Great Job!" is displayed. Only two acceptable, reproducible spirometry maneuvers need be conducted in Frontline mode.

## 6.2 Test type 'FVL (Inspiration and Expiration)'

With this test mode, a deep inhalation follows the exhalation maneuver directly. Proceed in the same way as with the above-described FVC test. However, instruct the patient not to remove the spirette™ from his or her mouth after exhaling, but to follow up with a deep, maximum inhalation. Three acceptable tests should be conducted with this test as well.

## 6.3 Test type 'Slow VC'

Slow spirometry serves to determine the vital capacity and the lung volumes (see Chapter 6.8). You can repeat the maneuver several times. The best test is saved. Proceed as follows:

- Insert a spirette™ into the instrument. Ensure that the arrow on the spirette™ lines up with the arrow on the instrument.
- Press ENTER when the patient is ready. You will now hear the sensor buzz.
- The instrument now prompts you to avoid flow in the spirette™ while it is setting the baseline. It is advisable to block off the spirette™ at one end to ensure the baseline is set precisely. An audible signal sounds when the zero point is set.
- The patient must now insert the spirette™ into his or her mouth and breathe at rest (about 3-5 breaths) until you hear an audible signal.
- The patient must then take a deep inspiration followed by a maximum exhalation.
- The instrument stops automatically at the end of the maneuver.

If you are only interested in the vital capacity without determination of the other volumes the VC maneuver can also be performed without waiting for the acoustic signal.

At the end of the SVC test you can decide to immediately add an FVC test. If you do so the parameter FEV1/VC (Tiffeneau) is also shown on the printed report of the FVC test.

## 6.4 Test type 'MW'

- Insert a spirette™ into the instrument. Ensure that the arrow on the spirette™ is lined up with the arrow on the instrument.
- Press ENTER when the patient is ready. You will now hear the sensor buzz.
- The instrument now prompts you to avoid flow in the spirette™ while it is setting the baseline. It is advisable to block off the spirette™ at one end to ensure precise setting of the baseline. An audible signal sounds when the baseline is set.
- The patient must then insert the spirette™ into his or her mouth and must fully inhale and exhale for an uninterrupted period of at least 12 seconds.

## 6.5 Test type 'OSHA Cotton Dust'

This test type is available only if 'North America' has been selected as the region (see Chapter 4.3). The procedure for conducting the test is equivalent to that for test type 'FVC (Expiration)'. The quality criteria applied and the printout correspond to the specifications of the American Occupational Health and Safety Administration (OSHA) [2].

## 6.6 Test type 'Disability'

This test type is available only if 'North America' has been selected as the region (see Chapter 4.3). The procedure for conducting the test is equivalent to that for Test type 'FVC (Expiration)'. A calibration check is required before conducting the actual spirometry test. The quality criteria applied and the printout correspond to the specifications of the American Social Security Administration Disability (SSD) Guidelines [3].

## 6.7 'Post' test

The Post-Test is usually performed to determine the response to bronchodilating asthma medication. This is done by treating the patient with a bronchodilator after having performed a FVC or FVL test. Approximately 10 to 20 minutes after the medication (when bronchodilator shows effect) a second FVC or FVL test ("post-Test") is performed. The results of the pre-test and the post test are then compared on the result screen and on the test report. Post-Tests can only be added to previous tests on the same day.

## 6.8 Quality messages

The table below lists all quality messages that may be displayed at the end of a test:

Note	Note relates to	Suggestion for improvement
Don't hesitate	...quality of the last test	The patient should exhale all at once, and not in steps.
Breathe out faster	...quality of the last test	The patient must exhale explosively, as hard and quickly as possible
Breathe out longer	...quality of the last test	The patient has stopped exhaling too early. The patient must exhale still further and force as much air as possible out of his or her lungs.
Good effort, do next	...quality of the last test	Good try. Just one to two further good tries and the test is complete.
Breathe out harder (only in Frontline mode)	...reproducibility of the spirometry maneuvers. PEF not reproducible	The test differs greatly from the previous tests. The patient can breathe out even harder and achieve a higher peak flow.

Breathe more deeply	...reproducibility of the spirometry maneuvers. FVC or FEV1 not reproducible	The test differs greatly from the previous tests. The patient can inhale even more deeply and exhale even more air.
Test complete	...quality of the entire test	The test is complete. An adequate number of good spirometry maneuvers has been achieved.

## 6.9 Interpretation of results

The quality grades allow you to assess the reliability of the measurement result. Quality grades A to C indicate a reliable result. A quality grade of D or F indicates inadequate test quality. The result must then be interpreted with caution. You can either print out the report immediately with (ENTER) or click on the DATA button to view the result on the display as soon as you see message “Session Complete! Great Job!” after performing a test.

Measured values below the lower limit values are shown in red on the report and explicitly indicated with \* and a footnote. EasyOne-CS also offers an automatic interpretation aid. Please refer to Chapter 0 for further information on this interpretation function. It is possible to deactivate both the “Automated QC” and “Syst. Interpret” functions.

## 6.10 Definition of the parameters

Below is a list of available parameters, with their definitions:

FVC	Forced Vital Capacity (expiratory)
FIVC	Forced Vital Capacity (inspiratory)
FEV1	Forced Expiratory Volume (1 sec)
FEV6	Forced Expiratory Volume (6 sec)
FEV1/FVC	Ratio of FEV1 to FVC
FEV1/VC	Ratio of FEV1 to VC taken from SVC test
FEV1/FEV6	Ratio of FEV1 to FEV6
MEF 25	Mid Expiratory Flow at 75% of exp. Vital Capacity
MEF 50	Mid Expiratory Flow at 50% of exp. Vital Capacity
MEF 75	Mid Expiratory Flow at 25% of exp. Vital Capacity
MEF 25-75	Mid Expiratory Flow between 75% and 25% of exp. Vital Capacity
PEF	Peak Expiratory Flow
PIF	Peak Inspiratory Flow
FET	Forced Expiratory Time
PRE/POST %	Percentage variation of measured values before and after bronchial spasmolysis
LLN	Lower Limit of Normal
BEV	Back Extrapolated Volume
VT	Tidal Volume
ERV	Expiratory Reserve Volume
IRV	Inspiratory Reserve Volume
VC or VCmax	Maximum Vital Capacity
VCex	Expiratory Vital Capacity
VCin	Inspiratory Vital Capacity
IC	Inspiratory Capacity
MVV	Maximum Voluntary Ventilation
Lung Age	See Chapter 16, [8]


## 7 EasyWare PC software

The EasyWare PC software can also be used independently of CardioSoft™. The software is used to view and print tests conducted with EasyOne.

**Please note:** Tests conducted only with EasyOne-CS and EasyWare (without CardioSoft™) cannot be imported to CardioSoft™ subsequently.

### 7.1 EasyWare's main window

The illustration below shows EasyWare's main window.

From top to bottom, you can see menus, toolbar, list of all records from the PC database and status bar.

The list can be sorted in ascending or descending sequence by record number, type, ID, name, age or test date by clicking on the corresponding column header.

The following information is available for each record, i.e. for each test conducted:

<b>Symbol</b>	Identification of the status of the record (see below)
<b>No.</b>	Internal record number. This is also displayed in EasyOne when selecting the test.
<b>Type</b>	Type of test conducted
<b>ID</b>	Patient identification number (if entered)
<b>Name</b>	Name of the patient (if entered)
<b>Age</b>	Age of the patient
<b>Test date</b>	Date and time of the test

The symbols below serve to identify the record status.

-  A red exclamation mark indicates a record that has not yet been synchronized with the EasyWare PC database, i.e. this record is currently saved only in EasyOne.
-  A blue database symbol identifies a record that is available in the EasyWare PC database.
-  A green database symbol indicates a record that has been generated via CardioSoft™ and is available in the EasyWare PC database. Certain operations, for example editing patient data, are not possible with this type of record.
-  A red database symbol indicates a record that exists in the PC database but that has been updated or supplemented in EasyOne. One example is adding a Post test. This record is synchronized during the next synchronization operation.
-  An amber sheet indicates a patient record that was generated in EasyWare. This record contains only patient data and no measurement data. This record can be used in EasyOne for future measurements.
-  A grey sheet indicates a patient record that was generated via CardioSoft™. This record contains only patient data and no measurement data.

The status bar, i.e. the lower area of the EasyWare main window, shows the serial number of the connected instrument and the current connection status (for example 'Not connected' or 'Device connected' etc.).

### Synchronization range

Records are within the synchronization range if their date is more recent or the same as the date of the most recent record in the database.

### Note

Note that the first synchronization operation for an EasyOne may take a relatively long time since, under certain circumstances, many records may have already been saved in EasyOne.

## 7.2 PC database


A Microsoft Access-compatible database saves all data of an EasyOne record including curve data. The default name of the database file is EasyWare.mdb. The path indicating where the file was saved can be specified in the program properties (File menu / Properties).

The EasyOne database is normally synchronized automatically with the PC database. All records with a red exclamation mark or a red database symbol are then saved in the database. It takes approx. 0.2 seconds to transfer one record.

After synchronization, the list of records is refreshed and only 'blue database symbols' or 'amber patient record symbols' are displayed.

## 7.3 Displaying patient data

Patient data can be displayed or edited using the Patient Data item in the Edit menu or using the  symbol in the toolbar. The following dialog box is displayed:

The 'Edit Patient Data' dialog box contains the following fields and controls:

- Name: Text box containing 'Mustermann, Manfred'
- ID: Text box containing '1'
- Birth Date: Date picker showing '1/ 2/1966'
- Height [ft/in]: Two text boxes containing '6' and '2'
- Weight [lbs]: Text box containing '63'
- Gender: Dropdown menu showing 'Male'
- Asthma: Dropdown menu showing 'No'
- Ethnic: Dropdown menu showing 'CAUCASIAN'
- Smoker: Dropdown menu showing 'Never Smoker'
- Buttons: 'OK' and 'Cancel'

If the EasyOne is connected to the PC, patient data can also be edited on the instrument. The data in EasyOne is updated automatically.

If using CardioSoft™, patient data is transferred automatically to EasyWare. Patient data generated in this way can neither be changed nor deleted in EasyWare or EasyOne-CS.

## 7.4 Entry of patient data

Choose New Patient in the Edit menu or use the  symbol in the toolbar to create a new patient record. Enter the patient data and confirm the dialog box. A new patient record is now generated and is displayed at the top of the list.

If the 'Technician ID' option has been activated in EasyOne's Device Configuration settings (see Chapter 8.3 'System Configuration'), a corresponding ID can be entered.

Note that this function is available only if the EasyOne is connected to the PC.

Also note that, if using CardioSoft™, patient data is entered in CardioSoft™ and not in EasyWare.

## 7.5 Deleting records


Records can be removed from the database using Delete in the Edit menu. Selecting several records at once (by using the Shift or Ctrl keys) also allows you to delete a range of records.

Note that records cannot be deleted within the synchronization range (see above). These records are reloaded from the instrument during the next synchronization operation.

Records that have been created in conjunction with CardioSoft™ cannot be deleted.

## 7.6 Test result preview

Test data can be assessed by selecting a record and double-clicking on it. Alternatively, you can also use Test Results in the View menu or the corresponding icon in the toolbar. The following window is opened:

The Test Result Preview window has the following content:

- Measurement curves (depending on type of measurement), Pre and/or Post test with predicted points.
- Patient and test information.
- Test results, Pre and/or Post Test
- Test quality and system interpretation.

You can skip to the next or previous record using the arrow keys. You can also do this on the keyboard using the PgUp or PgDn keys.

More information can be shown using the 'Advanced' selection button. Advanced parameters are displayed and you have the option of adding comments and editing as follows:

- Entry of a reviewer name is necessary for making changes.
- A comment can be assigned to the test.
- The order of the tests and the acceptance can be changed by clicking on a test (for example, click on 'Trial1'). This is possible only if the test was conducted with setting '3 best curves' (see Chapter 8.2 'Test Settings').

Tests performed using CardioSoft™ can be changed only if the test is displayed via CardioSoft™ (see Chapter 5.2.3).

The EasyWare Test Result Preview can be copied easily to other applications (Microsoft Word for instance). Perform the following steps to do this:

1. Launch the program to which the data is to be copied (Microsoft Word for instance).
2. Launch EasyWare.
3. Choose the test to be copied and double-click on the test or choose Test Results in EasyWare's View menu.
4. Choose what data is to be displayed (Baseline, Post or both)

5. Press and hold key Alt. Now also press the PrtSc (PrintScreen) key. This copies the contents of the active window to the clipboard.
6. Switch to the application to which the data is to be copied.
7. Press Ctrl-V or choose Insert in the Edit menu. This copies the contents of the clipboard to the application.


## 7.7 Online test

The procedure for conducting online tests is described in Chapter 5.4. The procedure for conducting tests with CardioSoft™ is described in detail in Chapter 5.2.

## 7.8 Pediatric animation

If using the EasyOne Screen-Connector, you can choose pediatric animation for online FVC tests. The following illustrations show the animation during the test and at the end of a good test.

A bursting balloon signifies a good test

Note that the pediatric animation is not available for test types FVL, SVC test or MVV.


## 7.9 Trend display

Clicking on Trend in the View menu allows you to display the trend of test data. The tool can be used in the following way:

- Sort the list by patient name or patient ID by clicking the corresponding header.
- Select the tests to be used for the trend display. Use key Shift or Ctrl to select more than one test.
- Select Trend in the View menu.

A trend display is now shown in the window. The number and type of parameters shown can be set at the left-hand side. The trend display can be sorted by date/time or by record number.

## 7.10 Searching for records

You can search for records by name or ID using Search in the View menu or by clicking on the corresponding symbol in the toolbar. Entering part of a name or number is also possible (for example “Thomas” is also found if you enter “Thom”). The search always begins at the current position.

You can continue the search by clicking on Find Next in the View menu.

## 7.11 Printing a test

Choose Print Preview in the File menu (or click on the corresponding symbol) to display a preview. Choose Print in the File menu (or click on the corresponding symbol) to print the test. The type of print-out depends on the EasyOne Device Configuration settings (see Chapter 8.4).

## 7.12 Exporting records

EasyWare supports the following data export formats:

Text export (CSV): All Records	All records of the current instrument are exported as text files. CSV (Comma Separated Values) is used as the data format. See Appendix B for details. The file names are generated based on the instrument serial number and the record numbers.
Text export (CSV): New Records	New records of the current instrument are exported as text files. CSV is used as the data format. See Appendix B for details. The file names are generated based on the instrument serial number and the record numbers.
Text export (CSV): All Devices	All records of all existing instruments are exported as a CSV text file. A test is shown on one single line in the text file. This line contains the patient data and selected data of the test.
GDT export	A GDT file is generated on the basis of the selected record. GDT files are used for exchange with patient information systems.
XML export (XML)	The selected records are exported as XML files. The file names are generated based on the instrument serial number and the record numbers.
ATS/ERS export (CSV)	The selected records are exported to a CSV file. The data format corresponds to the ATS/ERS format (ATS/ERS Task Force: Standardization of Lung Function Testing. Eur Respir J 2005; 26:319-338).

The files are exported to the folders defined in Properties in the File menu.

### 7.13 Importing records

EasyWare allows import of XML, p7m (encrypted e-mail) or eml (unencrypted e-mail) files.

### 7.14 Selecting the EasyOne instrument for the database view

As mentioned above, several EasyOne instruments can be synchronized with one PC database. The individual instruments are distinguished in the database by their serial number. In Offline mode, i.e. if no instrument is connected to the software, the required instrument serial number can be selected manually. Clicking on Select Device Database in the View menu selects the instrument. All instrument serial numbers present in the database are shown in a list. Now choose the required serial number.

### 7.15 Device status


The device status (Device Status in the View menu) provides information on serial number, instrument software version, instrument boot software version, hardware version and device type of EasyOne.

### 7.16 Device Configuration

The entire Device Configuration (under Device Configuration in the EasyOne main menu) can be changed easily with EasyWare. The possible settings are identical to the settings in EasyOne. For instance, the header of the report can be changed easily. The clock in EasyOne can be synchronized with the PC clock by clicking on a button on the 'System Configuration' page.

### 7.17 Program Preferences

The following window is displayed if you click on Properties in the File menu.

The following properties can be changed:

- Type of connection to PC. The 'Screen-Connector' is selected here automatically for EasyOne-CS.
- File name and path of the EasyWare database. The default setting is EasyWare.mdb and the path where EasyWare is installed.
- Path for text file export (default setting the same as for database).
- Path for the GDT interface. These values are automatically set correctly for CardioSoft™.

### 7.18 Language Selection

The language of the EasyWare application can be selected by clicking on Select Language... in the File menu. Note that this selected language is independent of the language selected in EasyOne.

## 8 Device Configuration

The device configuration can be selected in the EasyWare software (Device Configuration in the View menu) or directly in EasyOne. If you wish to change the device configuration directly in EasyOne, choose the DEVICE CONFIGURATION option in the main menu. You will now see the Configuration menu. The tables below provide an overview of the setting options offered by EasyOne-CS. Choose the corresponding option as required.

### 8.1 Operating modes 'Diagnostic', 'Frontline' and 'NLHEP.'

The EasyOne-CS spirometer has three operating modes to choose from:

- **Diagnostic** mode offers you all available tests and setting options. Spirometry tests can be conducted in accordance with the guidelines of the ATS (American Thoracic Society) and ERS (European Respiratory Society).
- **NLHEP** mode allows you to conduct spirometry tests in accordance with the greatly simplified recommendations of the 'National Lung Health Education Program' [4].
- **Frontline mode** allows you to conduct simple expiration tests quickly. The NLHEP recommendations have been broadened slightly for this.

The differences between the operating modes are described in the table below.

	Diagnostic	Frontline	NLHEP
Test types	FVC (Expiration), F/V Loop (Inspiration and Expiration), slow VC, MVV, Pre-Post measurement	FVC (Expiration), Pre-Post measurement	FVC (Expiration), Pre-Post measurement
Parameters	FEV1, FVC, FEV1/FVC, FEV6, FEV1/FEV6, MEF25-75, MEF25, MEF50, MEF75, PEF, FET, FVC, PIF, IVC, IRV, ERV, FEV1/VC, MVV, Pre-Post % variation, QC level	FEV1, FEV6, FEV1/FEV6, FVC, FEV1/FVC, PEF, Pre-Post % variation, QC level	FEV1, FEV6, FEV1/FEV6, pre-post % variation, QC level
Quality control	Requires 3 acceptable, reproducible spirometry maneuvers. Details in Chapter 9	Requires 2 acceptable, reproducible spirometry maneuvers. Details in Chapter 9	Requires 2 acceptable, reproducible spirometry maneuvers. Details in Chapter 9
Automated QC	Quality control can also be defined differently manually	Automated QC is always active	Automated QC is always active
Test display and storage	Display and storage of best or 3 best tests including curves	Only display and storage of the best test	Only display and storage of the best test
Report configuration	User can select curve display and size	Predefined report with small FV and VT curves	Predefined report with small FV and VT curves

The default setting of the EasyOne-CS Diagnostic Spirometer is Diagnostic mode. Please read Chapter 8.3 for information on how to change the operating mode.

### 8.2 Test Settings

Relates to	Option	Default setting	Description
Predicted	(1)	ERS / Zapletal North Am.: NHANES III	You can select your desired predicted values from the predicted publications listed.
Additional pediatrics	(1)	None	You have the option of selecting different predicted values for children than those for adults.

Relates to	Option	Default setting	Description
Best value selection	Best value, best trial	Best trial North Am.: Best value	In "best value" setting, the relevant, best value from different tests is selected. "Best trial" selects the test which has provided the best results (see Chapter 9.3).
System interpretation	NLHEP, GOLD/Hardie, OFF	NLHEP	Automatic interpretation (see Chapter 0) can be activated or deactivated here.
Lung Age	yes, no	No North Am.: Yes	if set to "yes", the lung age is displayed on the result screen and printed on the report. Lung Age is only shown if the patient is a smoker.
Automated Test QC	Yes, no	Yes	The automated test QC (see also Chapter 9) is activated and deactivated here.
FVC selection	FVC, FEV6	FVC	FEV6 indicates the exhaled volume after 6 seconds. When set to FEV6, EasyOne stops the measurement after 6 seconds. MEF25, MEF50, MEF75 and MEF25-75 are not reported in that setting. When set to FVC, EasyOne continues the measurement until end of test criteria are met.
PEF unit	L/s, L/min, OFF	L/s North Am.: L/min	Peak flow can be specified in liters per minute or in liters per second. OFF: PEF is not shown.
African ethnic correction	75 % - 110 %	88 %	The predicted normal is corrected by this additional factor if no specific calculation is defined for this ethnic group in the selected predicted normal publication.
Asian ethnic correction	75 % - 110 %	100 %	The predicted normal is corrected by this additional factor if no separate calculation for this ethnic group is defined in the selected predicted normal publication.
Latin American ethnic correction	75 % - 110 %	100 %	The predicted normal is corrected by this additional factor if no specific calculation for this ethnic group is defined in the selected predicted normal publication.
Other ethnic correction	75 % - 110 %	100 %	The predicted normal is corrected by this additional factor if no specific calculation for this ethnic group is defined in the selected predicted normal publication.
Curve storage	3 best curves, best curve	3 best curves	EasyOne can save the 3 best curves of a test. This is necessary if you want to print the 3 best curves or if you want to export the curve data of the three best trials. Please note that saving 3 curves uses substantially more memory. EasyOne can then only save about 250 tests before it starts to overwrite the oldest test.

(1) The predicted normals available depend on the region selected.

### 8.3 System Configuration

Relates to	Option	Default Setting	Description
Time format	24 hours, am/pm	24 hours North Am.: am/pm	Sets the time format for 12 or 24 hour.
Date format	DD.MM.YY, DD/MM/YY, MM/DD/YY	DD.MM.YY North Am.: MM/DD/YY	Sets the data format.

Relates to	Option	Default Setting	Description
Current date			Please enter the correct date at this point and confirm with ENTER.
Current time			Please enter the correct time at this point and confirm with ENTER.
Alphanumeric ID	Yes, no	No	If the ID you use also consists of letters, please set to "Yes".
Technician ID	Yes, no	No	If you want the technician ID to be saved as well and listed on the report, please choose "Yes".
Syringe volume	1.0l, 1.5l, ...7.0l	3.0l	Choose the volume of your syringe if you wish to conduct a calibration check.
Height unit	m/cm, ft/inch	m/cm North Am.: ft/inch	Choose how the unit will indicate height and altitude.
Weight unit	kg, lbs	Kg North Am.: lbs	Choose how the unit will indicate weight.
Age, date of birth	Age, birth	Birth North Am.: age	If you use a database, consider entering the date of birth so that the age can be calculated correctly at a later point.
Contrast			Changes the display contrast.
Lang.	German, English, others	English	Choose the desired language.
Altitude above sea level	0m, ...4000m	0m	Set the altitude above sea level of your location.
Op. mode	Diagnostic, Front-line, NLHEP	Diagnostic	see Chapter 8.1
Temp. unit	°C, °F	°C North Am.: °F	Determines how temperature is reported.
Rel. humidity	0...100%	40%	Enter the average relative humidity at your location.

## 8.4 Print Settings

Relates to	Options	Default Setting	Description
Printer type	via PC	via PC	Set to 'via PC' by default if using CardioSoft™.
Result data	3 best values, best values	3 best values	You have the choice of printing out only the best test or the 3 best tests on the report.
Number of Curves	3 best curves best curve	best curve	Choose if you want to print the 3 best curves of the tests or only the best curve. It is only possible to print the 3 best curves if the 3 best curves were saved (see test settings "Save curve data").
Graph Types	FV&VT small, FV large, VT large, FV&VT large	FV&VT small	Choose what curves you wish to have on the report.
Header 1-4	Optional entry	Blank	You can enter the name and address of the institution or other information in 4 lines of 40 characters each.

## 9 Quality messages and quality grades

### 9.1 Quality messages

The quality messages serve to provide assistance when conducting the measurement. After a test, they provide information as to whether the test was acceptable or what action can be taken to achieve a better result.

Message	Criterion	Recommended action
Don't hesitate	Back-extrapolated volume greater than 150 ml or 5% whichever is greater	The patient must exhale all air at once and not exhale in short bursts.
Breathe out faster	Time until peak flow greater than 120 ms	The patient must exhale more explosively and as firmly and quickly as possible.
Good effort, do next	Test meets above criteria	Good test. Only one to two more good tests and the test is complete.
Breathe out harder (only in Frontline mode)	Peak flow not reproducible. Difference with respect to best test greater than 1.0 l/s	The test differs greatly from previous best test. The patient can blow even more firmly and achieve a higher peak flow.
Breathe more deeply	FEV1 or FVC* not reproducible. Difference with respect to best test greater than 150 ml	The test differs greatly from previous tests. The patient can inhale even more deeply and exhale even more air.
Test complete	Three acceptable tests, FEV1 and FVC* within 200 ml / 250 ml (after 5 trials)	The test is complete. An adequate number of good tests are available.

\* when using FEV6 instead of FVC, FEV6 is also used for the determination of the quality message

### 9.2 Quality Grades

The quality grades allow you to assess the reliability of the measurement result. Quality grades A to C indicate a reliable result. Quality grades of D or F indicate inadequate test quality. The result must then be interpreted with caution.

The quality ratings can be activated or deactivated under “Configuration”. See also Chapter 8. The table below defines the criteria for the quality grades:

Rating	Criteria in Diagnostic mode	Criteria in Frontline mode
A	At least 3 acceptable tests AND the difference between the best two FEV1 and FVC values is equal to or less than 150ml	At least 2 acceptable tests AND the difference between the two FEV1 and FEV6 values is equal to or less than 100ml (or FVC is equal or less 150ml)
B	At least 3 acceptable tests AND the difference between the best two FEV1 and FVC values is equal to or less than 200ml	At least 2 acceptable tests AND the difference between the two FEV1 and FEV6 values is equal to or less than 150 ml (or FVC is equal or less 200 ml)
C	At least 2 acceptable tests AND the difference between the best two FEV1 and FVC values is equal to or less than 250ml	At least 2 acceptable tests AND the difference between the two FEV1 and FEV6 values is equal to or less than 200 ml (or FVC is equal or less 250 ml)
D	At least 2 acceptable trials but the results are not reproducible. Quality message "Result not reproducible" OR only one acceptable trial. Quality message: "Only one acceptable trial"	At least 2 acceptable trials but the results are not reproducible. Quality message "Result not reproducible" OR only one acceptable test. Quality message "Only one acceptable trial"
F	No acceptable test available	No acceptable test available

If the “Automated QC” function is activated, the instrument automatically determines test acceptability. The acceptable tests are used ahead of non-acceptable tests for best test, system interpretation and the Pre/Post comparison.

In the Diagnostic mode, the “Automated QC” function can be deactivated (see Chapter 8 'Device Configuration'). In this mode, the EasyOne’s determination of acceptability can be overridden by choosing either KEEP or REJECT after the test. KEEP marks the test acceptable and REJECT marks the test unacceptable.

### 9.3 Best Test Selection


In the system configuration the best value selection can be set to “Best Trial” or “Best Value”. The two settings are defined as follows:

- Best trial:** EasyOne-CS selects the best trial by searching the trial with the largest sum of FVC and FEV1 (this is suggested by ATS and ERS). An acceptable test is always used ahead of an unacceptable test.
- Best value:** The “Best” column shows the largest FVC (or FEV6) and the largest FEV1 from all acceptable tests (unless all tests are unacceptable). All other parameters are taken from the best trial (again defined by the largest sum of FEV1 and FVC).

## 10 System Interpretation

The diagrams below describe how the interpretation is determined (see references [4], [12],[13]). Automatic interpretation (Setting: NLHEP or optional GOLD/Hardie) can be activated or deactivated under “Configuration” (see Chapter 8).


### 10.1 NLHEP System Interpretation

- Notes:
1. LLN = Lower Limit of Normal
  2. Where FVC is indicated in chart FEV6 may be substituted, if used.
  3. Where there is no lower limit of normal (LLN) defined within the selected predicted normal, the value used for LLN is calculated as Predicted Value – 1.645 x SEE (standard error of the estimate). If SEE is not defined LLN of FEV1/FVC is set to 90% of Predicted Value, LLN of FEV1 is set to 80% of Predicted Value; LLN of FVC is set to 80% of Predicted Value.
  4. If the quality grade is D and the results are within normal limits, the interpretation states “normal, but the reported values should not be used for comparisons with previous or subsequent tests.”

## 10.2 GOLD/Hardie System Interpretation

## 11 Predicted Values

EasyOne-CS offers a number of published predicted value tables allowing comparison of the measurement results. In order to compute the predicted values, it is necessary to enter the sex, age and height and, in many cases, the ethnic group and the weight of the patient. See also Chapter 8 on selection of the predicted values.

Where there is no lower limit of normal (LLN) defined within the selected predicted normal, the value used for LLN is calculated as Predicted Value – 1.645 x SEE (standard error of the estimate). If SEE is not defined LLN of relational parameters, e.g. FEV1/FVC are set to 90% of Predicted Value, LLN of all other parameters are set to 80% of Predicted Value.

If the patient data lies outside of the range defined in the publication (Age, Height) EasyOne-CS uses extrapolated values. The report explicitly points out that the predicted values are extrapolated values and that, consequently, particular caution must be taken when interpreting the results.

### 11.1 Predicted Values for Adults

- [1] Morris, James F., Koski, Arthur, Lavon Johnson, "Spirometric Standards for Healthy Non-Smoking Adults", *Am Rev Resp Dis*, Vol. 10-3, 1971. p.57-67
- [2] Morris, James F., Koski, Arthur, Breese, John, 'Normal Values and Evaluation of Forced Expiratory Flow', *Am Rev Resp Dis*, Vol. 111, 1975, p.755-761.
- [3] Cherniak, R.M. and Raber M.B. 'Normal Standards for Ventilatory Function using an Automated Wedge Spirometer', *Am Rev Resp Dis*, Vol. 106, 1972, p.38-46.
- [4] Knudson, Ronald J., Ronald Slatin, Michael Lebowitz, Benjamin Burrows, 'The Maximal Expiratory flow-Volume Curve', *Am Rev Resp Dis*, Vol. 113, 1976, p. 587-600.
- [5] Knudson, Ronald J., Michael Lebowitz, Holberg Catherine J., Benjamin Burrows, 'Changes in the Normal Maximal Expiratory Flow-Volume Curve with Aging', *Am Rev Resp Dis*, Vol. 127, 1983, p. 725-734
- [6] Crapo, Robert O., Gardner Reed M., Morris Alan H., 'Reference Spirometric Values Using Techniques and Equipment that Meets ATS Recommendations', *Am Rev Resp Dis*, Vol. 123, 1981, p. 659-674
- [7] Forche G., Harnoncourt K., Stadlober E.. 'Neue spirometrische Bezugswerte für Kinder, Jugendliche und Erwachsene', *Öst. Ärztztg.* 43/15/16 (1988) 40
- [8] P.H. Quanjer. Lung Volumes and Forced Ventilatory Flows. *Eur Respir J*, Vol 6, Suppl 16, p.5-40, 1993
- [9] NHANES III: Hankinson, Odencrantz, Fedan, 'Spirometric Reference Values from a Sample of the General U.S. Population,' *Am J Respir Crit Care Med*, Vol. 159, 1999, p 179-187.

### 11.2 Predicted Values for Children

- [1] Zapletal, T. Paul, M. Samanek. "Die Bedeutung heutiger Methoden der Lungenfunktionsdiagnostik zur Feststellung einer Obstruktion der Atemwege bei Kindern und Jugendlichen." *Z. Erkrank. Atm.-Org.*, Volume 149, 343-371, 1977
- [2] Hsu KHK, Bartholomew PH, Thompson V, Hseih GSJ, "Ventilatory Functions of Normal Children & Young Adults – Mexican-American, White, Black. I. Spirometry", *J. Pediatrics* 1979, 95: p.14-23
- [3] NHANES III: Hankinson, Odencrantz, Fedan, 'Spirometric Reference Values from a Sample of the General U.S. Population', *Am J Respir Crit Care Med*, Vol. 159, 1999, p 179-187.
- [4] Dockery, D.W. et al., Distribution of Forced Vital Capacity and Forced Expiratory Volume in One Second in Children 6 to 11 Years of Age, *Am Rev Resp Dis*, Vol. 128, p. 405-412, 1983
- [5] Polgar, Promadhat, Pulmonary Function Testing in Children: Techniques and Standards. W.B. Saunders Co., Philadelphia, 1971

## 12 Hygiene and instrument maintenance

When used correctly the EasyOne-CS requires minimal maintenance and service.

When you use the spirette™ respiratory tube, you do not need to clean the instrument. Instead of cleaning, you simply change the respiratory tube. In order to ensure absolute hygiene, we recommend that the spirette™ be used only once per patient.

**Caution:** Always exchange the spirette™ if you suspect the risk of infection. This is the only way of absolutely preventing transmission of diseases.

Use a damp cloth to clean the spirometer housing and the base unit. Use a soft cloth and alcohol (e.g. isopropyl alcohol) for particularly thorough cleaning.

**Caution:** Avoid fluid penetrating the spirette™ holder or the inside of the instrument when cleaning the spirometer.

To change the batteries please proceed as follows: Open the battery lid on the back side of the instrument. Remove the old batteries. Insert two new alkaline batteries (Type AA, 1.5V). Put the battery lid back on. Ensure that the battery compartment cover is firmly closed since patient safety may depend on it.

In the event of instrument malfunction, please contact your EasyOne-CS dealer or call the After-Sales Servicing Department of GE Medical Systems Information Technologies (<http://www.gehealthcare.com>).

Please proceed as follows if you wish to check whether your instrument is functioning correctly:

1. Check calibration. Please read Chapter 13.
2. Conduct a spirometry test on yourself.
3. Ensure that the results are plausible and that you can print out the report as required.

Consult your EasyOne dealer if you encounter problems with one of these points.


## 13 Calibration Check

Calibration of the instrument can be checked with a syringe using the calibration check function. The American Thoracic Society (ATS) recommends that calibration be checked every day. The EasyOne-CS spirometer requires no calibration, even if used frequently.

To perform a calibration check, you need the optional ndd calibration adapter and an optional calibration syringe (order number 2030431-010), in addition to the spirometer and a spirette™. Ensure that the correct syringe volume is entered in the instrument's configuration setting (see also Chapter 8).

Proceed as follows:

- Choose item "Check calibration" in the main menu.
- Connect the spirometer as shown below using the calibration adapter and the syringe. Ensure that the piston is fully inserted and at the stop position.

- Now press ENTER
- Wait until the baseline has been set and you hear an audible signal.
- Now execute one full inspiratory pump stroke followed by one full expiratory pump stroke at moderate speed.
- After you perform the maneuver, you will see the text "Accuracy confirmed" at the top of the display and, beneath it, the percentage deviation and the average flow velocity of the pump stroke.
- You can repeat the test, print the result or quit the program. The calibration test remains stored and can also be viewed or printed out later.

If you do not reach  $\pm 3\%$  accuracy, please follow the troubleshooting instructions in Chapter 14. Please contact the After-Sales Servicing Department at GE Medical Systems Information Technologies (<http://www.gehealthcare.com>) if you are unable to remedy the fault or error even by following this information.

## 14 Troubleshooting Tips

The table below provides a few tips for simple troubleshooting, should you have problems operating your spirometer.

Problem	Possible cause	Solution
EasyOne-CS cannot be switched on	Batteries are dead	Insert new batteries.
	Batteries are inserted wrong	Insert the batteries correctly (see Chapter 4.1.)
	The ON/OFF key has not been pressed for at least 2 seconds	Press the ON/OFF key for at least 2 seconds.
Three consecutive tones sound as a warning signal when the EasyOne-CS is switched on	The spirometer is defective	Consult your EasyOne dealer.
A self-test error message is displayed when the EasyOne-CS is switched on	The spirometer may be defective	Switch off the instrument and switch it back on again. Try again. If the error reoccurs, please contact your EasyOne dealer.
You are prompted to enter the date etc. each time you switch on	The internal battery of the EasyOne-CS is defective	Consult your EasyOne dealer.
The following message is displayed after a test is started: "Please insert spirette™ correctly"	The spirette™ is not correctly positioned	Ensure that the triangle on the spirometer is lined up with the triangle on the spirette™.
EasyOne-CS is out of the $\pm 3\%$ range in the calibration check	The spirette™ is not correctly positioned	Insert the spirette™ as described in Chapter 4.1.
	You have not used an ndd adapter	Use the ndd calibration adapter.
	There are leaks in the syringe connection	Check the connections.
	The specified syringe volume does not correspond to the actual syringe volume	Choose the correct syringe volume under "Configuration".
There is no curve on the printout	The color cartridge of your printer is empty	Replace the cartridge.
Only meaningless characters are printed out or the printer is not responding	The wrong printer type has been selected in Settings	Choose the right printer in the PRINT SETTINGS menu.
	The printer cable is not connected directly or defective	Switch off the spirometer and the printer. Check all plug connections.

Problem	Possible cause	Solution
	The printer is not switched on or not ready	Ensure that the printer is switched on and that it also has paper. Switch the printer off again and back on again.
	EasyOne-CS is not connected correctly to the Screen-Connector	Connect the Screen-Connector correctly to the EasyOne-CS.
An instrument self-test error #20 is displayed when you switch on the EasyOne-CS	A spirette™ has been inserted during the switch-on process OR it is not inserted correctly	Switch the instrument back on again with spirette™ inserted or without spirette™. If the error message is repeated, please contact your dealer.
An instrument self-test error #14 or #15 is displayed when you start a new test.	The spirette™ is not positioned correctly	Insert the spirette™ as described in Chapter 4.1 and try again. If the error message reoccurs, please contact you EasyOne dealer.
In instrument self-test error #25 is displayed when you switch on the EasyOne	The internal battery of the EasyOne-CS may be defective	Switch the instrument back on again. If the error message reoccurs, please contact your EasyOne dealer.
The instrument has been configured for the wrong region.	The wrong region has been selected when you switched on for the first time.	Perform a complete instrument reset: switch the instrument off and back on again and enter digit string 1-3-5. Press Enter to perform a full reset.

## 15 Bibliography

- [1] American Thoracic Society. Standardization of Spirometry: 1994 Update, Nov. 11, 1994. Am J Resp Crit Care Med 1995; 152:1107-1136.
- [2] Occupational Health and Safety Administration (OSHA), Pulmonary Function Standards for Cotton Dust, 29 CFR: 1910.1043 Appendix D.
- [3] Social Security Administration Disability (SSD) Guidelines, CFR404: Appendix 1 to Subpart P.
- [4] Ferguson GT, Enright PL, Buist AS, et al. Office spirometry for lung health assessment in adults: A consensus statement from the National Lung Health Education Program. Chest 2000; 117:1146-1161.
- [5] ATS Pulmonary Function Laboratory Management and Procedure Manual, American Thoracic Society, New York, NY 10019.
- [6] Enright PL, Hyatt RE. Office Spirometry. Lea & Febiger, Philadelphia, 1987.
- [7] Hyatt, RE, Scanlon PD, Nakamura M. Interpretation of Pulmonary Function Tests – A Practical Guide. Lippincott – Raven, Philadelphia, 1997.
- [8] American Thoracic Society. Lung Function Testing: Selection of Reference Values and Interpretative Strategies, Am Rev Respir Dis 1991; 144:1202-1218.
- [9] Morris JF, Temple W. Short Report: Spirometric „Lung Age“ Estimation for Motivating Smoking Cessation, Preventive Medicine 14. 655-662 (1985).
- [10] Polgar, Promadhat, Pulmonary Function Testing in Children: Techniques and Standards. W.B. Saunders Co., Philadelphia, 1971.
- [11] M.R. Miller, R. Crapo, J. Hankinson, et al. ATS/ERS Task Force: Standardization of Lung Function Testing. Numbers 1 to 5. Eur Respir J 2005; 26: 153-161, 319-338, 511-522, 720-735, 948-968.
- [12] Global Strategy For The Diagnosis, Management, And Prevention Of Chronic Obstructive Pulmonary Disease, Executive Summary, Updated 2003 (GOLD)
- [13] Hardie et.al., "Risk of over-diagnosis of COPD in asymptomatic elderly never-smokers" Eur Respir J 2002;20: 1117-1122

## 16 Accessories

2030434-001	EasyOne-CS Quick-Start Guide (GER, ENG, FRE, ITA, SPA, BRA)
2030434-002	EasyOne-CS User Manual (GER, ENG, FRE, ITA, SPA, BRA)
2030431-003	Screen-Connector (USB connection to PC)
2030431-005	CD EasyWare V2.9.4
2030431-006	spirette Standard 50 ea. box
2030431-007	spirette Standard 200 ea. box
2030431-008	Nose clips, 25 ea.
2030431-009	Nose clip pads, 100 ea.
2030431-010	Calibration pump with calibration adapter
73700008	Alkaline battery (2 ea. required)


**Please note:** To ensure patient safety and interference-free operation and to guarantee the specified measuring accuracy, we recommend only original equipment accessories as available through GE Medical Systems distribution. The user is responsible for the application of accessories from other manufacturers.

## 17 Appendix A: Specifications

### 17.1 EasyOne-CS Model 2001 Spirometer

Size:	83 x 158 x 43 mm
Weight:	242 g
Measuring accuracy:	Volume: $\pm 2\%$ or 0.050 l Flow: $\pm 2\%$ or 0.020 l/s, (apart from PEF) PEF: $\pm 5\%$ or 0.200 l/s MVV: $\pm 5\%$ or 5 l/min
Measuring range:	Volume: $\pm 12$ l Flow: $\pm 16$ l/s
Resistance:	Less than 0.3 cm H <sub>2</sub> O/l/s
Display:	64 x 160 graphic display
Data input:	14-key keypad
Data memory capacity:	Up to 700 tests (in the internal memory of the EasyOne-CS)
Test types:	FVC, FVL, slow VC, MVV, Pre/Post
Parameters:	FVC, MVV, FEV6, FEV1, FEV1/FVC, FEV1/VC, FEV1/FEV6, MEF25 (FEF75), MEF50 (FEF50), MEF75 (FEF25), MEF25%-75% (FEF25%-75%), PEF, FET, FIVC, PIF, IVC, ERV, IRV, Pre-Post % variation, lung age
Breathing tube:	Disposable spirette™ breathing tube
Measuring principle:	Ultrasonic transit time measurement
Adult normals:	NHANES-III, Knudson_83, Knudson_76, Crapo, Morris, ERS (ECCS/ECGS), Forche (Austria), Sapaldia (Switzerland), Roca (Spain), Berglund, Gulsvik, Hedenström, Gore, Chermiak (for MVV only)
Pediatric normals:	Dockery (Harvard), Hsu, Zapletal, Polgar, Hibbert
Power supply:	2 alkaline batteries, Type AA, 1.5 V
Power consumption:	Typically 0.6 W
Report format:	A4
Storage:	Temperature: -20 to 50 °C, relative humidity: 5 % to 95 %, Ambient pressure: 500 to 1060 hPa
Operating conditions:	Temperature: 0 to 40 °C, relative humidity: 5 % to 95 %, Ambient pressure: 500 to 1060 hPa
Approvals:	CE Declaration of Conformity C CSA US approval, CAN/CSA-C22.2 No. 601.1-M90, S1-94, CSA 601.1 Amendment 2:1998, UL Standard No. 2601.1. FDA 510(k) approval, K993921. EasyOne-CS complies with or excels the published targets of the European Respiratory Society (ERS), the American Thoracic Society (ATS) and the National Lung Health Education Program (NLHEP).
Device classification:	 Application section of Type BF In the case of operation with PC, safety class 2 in accordance with IEC/EN 60601-1 Device not suitable for use in flammable anesthetic gases in mixtures with O <sub>2</sub> or NO.

### 17.2 EasyOne-CS Model 2010 Screen-Connector

Size:	64 x 44 x 25 mm
Weight:	82 g
Power supply:	Via USB port
Power consumption:	Typically 0.15 W
Function:	Connects the EasyOne-CS spirometer to a PC
Connections:	Standard USB 1.1, Type A, connector for connection to PC
	 Pins: 1 = VBus, 2 = D-, 3 = D+, 4 = GND
Storage:	Temperature: -20 to 50 °C, relative humidity: 5 % to 95 %, Ambient pressure: 500 to 1060 hPa
Operating conditions:	Temperature: 0 to 40 °C, relative humidity: 5 % to 95 %, Ambient pressure: 500 to 1060 hPa
PC:	The PC must comply with the corresponding IEC Standard (e.g. IEC 60950-1). The user/operator must ensure that the requirements of IEC 60601-1-1 regarding safety of medical electrical systems are complied with.

## 18 Appendix B: Text file export format

The text files exported by EasyWare are named using the following conventions: SNsssss-rrrr.csv, whereby sssss stands for the unique serial number of the EasyOne and rrrr stands for the unique record number.

The following table describes record export in detail.

Parameter name	Suf- fix	Type	Description			
SerNr		int	EasyOne instrument serial number			
RecNum		int	Record number			
TypeOfTest		int	The following test types are defined: 2 = FVC                                3 = FVC including Post test 5 = FVC NLHEP                      6 = FVC NLHEP including Post test 8 = FVL                                9 = FVL including Post test 11 = MVV                               12 = Provocation 13 = Disability (USA)                14 = SVC test 15 = Calibration                     16 = Multi-flow calibration 20 = Disability (USA)                21 = Disability including Post test (USA) 23 = OSHA (USA)                    24 = OSHA including Post test (USA)			
NoOfTrials	_p	int	Number of tests, basic test			
FEV6Selected		boolean	Selection of FEV6 or FVC mode			
AutoQCOn		boolean	Selection of automatic quality check			
Storage Option		boolean	Selection of saving (only) best or 3 best curves			
Boolean1	_p	boolean	Not used			
Boolean2	_p	boolean	Not used			
Boolean3	_p	boolean	Not used			
Boolean4	_p	boolean	Not used			
Boolean5	_p	boolean	Not used			
QCGrade	_p	int	Quality grades: 0 = F (poorest quality)                1 = D 2 = D not reproducible                3 = C 4 = B                                        5 = A (best quality)			
Date		dd.mm.yy	Test date			
Word0	_p	word	Not used			
Word1	_p	word	Not used			
Word2	_p	word	Not used			
BtpsExp	_p	int	Expiration BTPS factor: 0.90 + BtpsExp/100			
BtpsIn	_p	int	Inspiration BTPS factor: 0.95 + BtpsIn/100			
PatientID		string	Identification number of the patient			
BirthDate		dd.mm.yy	Date of birth (if 'Age', then current date - age)			
Height		int	Size in cm			
Weight		int	Weight in kg*100			
Gender		int	Sex of the patient: 0 = male, 1 = female			
Smoker		int	Smoker status: 0 = yes, 1 = no, 2 = ex			
Asthma		int	Asthma status: 0 = no, 1 = possible, 2 = yes			
Ethnic		int	Ethnic group: 0 = African                                1 = Caucasian / European 2 = Latin American                      3 = Asian 4 = Other			
Name		string	Patient name			
TechnID		string	Identification number for technician / assistant			
Spare		int	Not used			
Accept	_p	3 int	Test quality: 0 = unacceptable, 1 = acceptable			
TrialTyp	_p	3 int	Not used			
TrialNo	_p	3 int	Test number			
Time	_p	3 hh:mm:ss	Time			
				<b>SVC test</b>	<b>MVV</b>	<b>Calib.</b>
FEV1	_p	3 float	Forc. exp. vol. in 1 sec. [l]	VT, tidal vol [l]		insp. Vol. [l]
FEV6	_p	3 float	Forc. exp. vol. in 6 sec. [l]	ERV [l]		exp. Vol. [l]
FVC	_p	3 float	Forc. vital capacity [l]	VCex [l]		Pump volume [l]
PEF	_p	3 float	Exp. peak flow [l/s]			
BEV	_p	3 float	Back-extrapolated volume [l]			
FET	p	3 float	Forc. exp. time [s]	SVC time [s]	MVV time [s]	

PEFT	_p	3 float	Exp. peak flow time [s]	Rf [1/min] / 2		
FEF25	_p	3 float	forc. exp. flow 25% [l/s]			
FEF50	_p	3 float	forc. exp. flow 50% [l/s]			
FEF75	_p	3 float	forc. exp. flow 75% [l/s]			
FEF2575	_p	3 float	forc. exp. flow 25-75% [l/s]			mean exp. flow
FEV3	_p	3 float	forc. exp. volume in 3 s [l]	IRV [l]	MVV [l/min]	
FIVC	_p	3 float	forc. insp. vital capacity [l]	VCin [l]		
VCmax	_p	3 float	max. vital capacity [l] (FVL only)			
EOTV	_p	3 float	Volume at end of test [l]			
PIF	_p	3 float	Insp. peak flow [l/s]	IC, insp.Cap.[l]		insp. Peak[l/s]
FIF25	_p	3 float	forc. insp. flow 25% [l/s]			
FIF50	_p	3 float	forc. insp. flow 50% [l/s]			
FIF75	_p	3 float	forc insp. flow 75% [l/s]			
Tzero	_p	3 float	Time datum for back-extrapolation [s]			
CuFVLen	_p	3 int	Length of flow/volume curve			
CuVTLen	_p	3 int	Length of flow/time curve			
CurveFV	_np	x int	Flow/volume curve			
CurveVT	_np	y int	Volume/time curve			

## Notes:

- In the case of test types ‘SVC test’, ‘MVV’ and ‘Calibration’, the content of certain parameters will differ from the conventional values. The corresponding values are noted at the end of the table.
- Curve data is exported in vertical format so that it can be imported to EXCEL easier.
- Boolean values: zero designates ‘False’ and not equal to zero designates ‘True’.
- Suffix p: this suffix is appended to a parameter name if the test is a Post test. Example: FEV1Best\_p designates the best FEV1 value in the Post test. parameters with suffix \_p exist only if a Post test has been conducted.
- Suffix \_np: this parameter suffix is used only for curve data. Up to three curves may be exported dependent on the Device Configuration: CurveFV\_1, CurveFV\_2, CurveFV\_3. The numbers signify 1 = best test, 2 = second best test, 3 = third best test. The additional suffix p is, in turn, used for Post tests.
- CurveFV contains precisely the number CuFVLen of values. The same applies analogously to CurveVT.
- Scaling of the FV curve: the flow values in the table have a resolution of 10 ml/s (a value of 10 means a flow of 100 ml/s). The volume spacing between two points is 30 ml. The first measured value at (0.0) is always omitted.
- Scaling of the VT curve: The volume values in the table have a resolution of 10 ml (a value of 100 corresponds to 1000 ml). The time spacing between two measurement points is 60 ms. The first point of the curve at (0.0) is always omitted.
- Parameters signs: inspiration positive, expiration negative.

## 19 Appendix C: Predicted Normals

This document describes in detail which predicted normals are implemented in the current version of EasyOne and the PC companion software EasyWare. The following table shows implementation details and the publications on which the reference values are based.

The first table has the following entries:

Reference:	Name of predicted normal reference used in EasyOne and EasyWare. The name is the study name, or the name of the authors of the publication.										
Publ. Year:	Year when the study was published.										
Abbrev.:	Abbreviation for the study. The abbreviation is used if parameters from other studies are copied. Example: For most studies 'Cherniak' (CH) values are used for the MVV parameter.										
Age Range:	Age range of the study. If range was extended beyond the study range the extended range is shown in parenthesis.										
Height Range:	Height range of the study. Extended range also listed in parenthesis, as in age range.										
Weight Range:	Weight range of the study. In most studies weight is not used.										
Built in:	The studies that are marked with an X are integrated into all EasyOne. Studies that are not marked can be loaded into EasyOne by ndd or a local dealer.										
Ethnic:	Ethnic group that is supported by the study.										
Parameter:	In the parameter list the following indicators are used: <table> <tr> <td>X</td><td>the parameter is available in the study</td></tr> <tr> <td>X (green)</td><td>the parameter is available in the study <u>and</u> it is used in EasyOne</td></tr> <tr> <td>X (orange)</td><td>the LLN of the parameter (see below) is computed according to the ATS recommendations: <math>LLN = Predicted - 1.645 * SEE</math> (Standard Error of Estimate)</td></tr> <tr> <td>FVC, IVC</td><td>the value of the parameter is copied from listed parameter of the same source</td></tr> <tr> <td>FEV1/FVC</td><td>FEV1/FVC% is computed from the predicted of FEV1 and FVC of the same source</td></tr> </table>	X	the parameter is available in the study	X (green)	the parameter is available in the study <u>and</u> it is used in EasyOne	X (orange)	the LLN of the parameter (see below) is computed according to the ATS recommendations: $LLN = Predicted - 1.645 * SEE$ (Standard Error of Estimate)	FVC, IVC	the value of the parameter is copied from listed parameter of the same source	FEV1/FVC	FEV1/FVC% is computed from the predicted of FEV1 and FVC of the same source
X	the parameter is available in the study										
X (green)	the parameter is available in the study <u>and</u> it is used in EasyOne										
X (orange)	the LLN of the parameter (see below) is computed according to the ATS recommendations: $LLN = Predicted - 1.645 * SEE$ (Standard Error of Estimate)										
FVC, IVC	the value of the parameter is copied from listed parameter of the same source										
FEV1/FVC	FEV1/FVC% is computed from the predicted of FEV1 and FVC of the same source										
Lower Limit of Normal of Parameter:	This table lists which lower limits of normal (LLN) are provided by the study. If the study does not have a formula for LLN, the LLN is set to 80% of the predicted for normal parameters (FVC, FEV1 etc.) and 90% of predicted for relative parameters like FEV1/FVC%.										

Additional remarks:

Lung Age: Lung Age is computed according to the following publication: Spirometric "Lung Age" Estimation for Motivating Smoking Cessation. J.F. Morris, W. Temple. Prev Med 14, 655-662 (1985).

Scandinavian References: These predicted normals are always combined with the Zapletal references for children.

								Ethnic					Parameter													Lower Limit of Normal of Parameter																								
	Reference	Publ. Year	Abbrev.	Age Range [years]	Height Range [cm]	Weight Range [kg]	Build In	Caucasian	African	Mexican	Asian	Other	FVC	FVC	VC	FEV1	FEV1/FVC%	FEV1/VC%	FEV3	FEV6	FEV3/FVC%	FEV1/FEV6%	FEF25	FEF50	FEF75	FEF25-75	PEF	MVV	Lung Age	FVC	FVC	VC	FEV1	FEV1/FVC%	FEV1/VC%	FEV3	FEV6	FEV3/FVC%	FEV1/FEV6%	FEF25	FEF50	FEF75	FEF25-75	PEF	MVV	Lung Age				
North America	NHANES III	1999	NH	8..80 (6..90)	110..195 (110..210)	unused	0	0	0				0	FVC	FVC	0	0			0		0				0		CH	0	0	FVC	FVC	0	0			0		0		0				0					
	Knudson_83	1983 (1976)	KN	6..85	107..183/196	unused	0	0					0	FVC	FVC	0	0		0	76					0	0	0	76	CH	0	0	FVC	FVC	0	0			-				0	0	0						
	Knudson_76	1976	KN	8..85	110..200 (110..210)	unused	0	0					0	FVC	FVC	0	0		0							0	0	CH	0	0	FVC	FVC	0	0																
	Crapo	1981	CR	15..91	146..195 (146..210)	unused	0	0					0	FVC	FVC	0	0		0		0						0		CH	0	0	FVC	FVC	0	0		0		0						0					
	Morris	1971 (1976)	MO	20..90	142..203 (142..210)	unused	0	0					0	FVC	FVC	0	76										0		CH	0	0	FVC	FVC	0	0									0						
	Hsu	1979	HS	7..20 (4..20)	110..195 (93..210)	unused	0	0	0				0	FVC	FVC	0	FEV1/FVC										0	0	CH	0	0	FVC	FVC	0									0							
	Dockery (Harvard)	1993	DO	6..18	115..185 (115..210)	unused	0	0	0				0	FVC	FVC	0	0										0		CH	0	0	FVC	FVC	0	0									0						
	Polgar	1971	PO	4..17	109..170 (90..210)	unused	0	0					0	FVC	FVC	0	0		0								0	0	0		0	FVC	FVC	0	0															
Cherniak	1972	CH	15..79	100..200	unused	0	0								0	0							0	0	0	0	0	0																						
Latin America	Pereira	1992	PE	20..78 (20..)	136..182 (90..220)	unused		0	0				0	FVC	FVC	0	0										0	0			0	FVC	FVC	0	0															
Europe	ERS (ECCS / EGKS)	1993	ER	18..70 (18..90)	145..195 (145..210)	unused	0	0					0	FVC	IVC	0	0						0	0	0	0	0	CH	0	0	FVC	IVC	0	0							0	0	0	0						
	Zapletal	1977	ZA	6..17 (4..17)	115..180 (93..210)	unused	0	0					VC	VC	0	0	0							0	0	0	0	0	CH		VC	VC	0	0	0						0	0	0	0						
	Forche (Austria)	1988	FO	7..76 (7..90)	110..200 (110..210)	unused	0	0					0	FVC	FVC	0	0						ER	ER	ER	ER	ER	0		0	FVC	FVC	0	0							ER	ER	ER	ER						
	Sapaldia (Swiss)	1996	SA	18..60 (18..90)	...	unused		0					0	FVC	FVC	0	0						0	0	0	0	0			0	FVC	FVC	0	0							0	0	0	0						
	Roca (Spain)	1982	BA	6..70	110-200 (110..210)	> 0	0	0					0	FVC	FVC	0	FEV1/FVC							0	0	0	0	CH																						
Europe Scandinavia	Hedenström + ZA	1985-1986	HE	20..70 (18..90)	150..195 (150..210)	female: 45..94	0	0					0	FVC	0	0	FEV1/VC	0					0	0	0		0		0	0	FVC	0	0	FEV1/VC	0						0	0	0							
	Gulsvik + ZA	1985	GU	20..70 (18..90)	150..190 (150..210)	unused	0	0					VC	VC	0	0	0									0		0		0	VC	VC	0	0	0								0							
	Berglund + ZA	1963	BE	20..70 (18..90)	154..191 (154..210)	unused	0	0					VC	VC	0	0	0											0		0	VC	VC	0	0	0															
Australia	Hibbert	1989	HI	8..19 (6..19)	120..190 (90..220)	unused							VC	VC	0	0							0	0	0	0				VC	VC	0	0									0	0	0	0					
	Gore, Crockett	1995	GO	18..78 (18..)	145..195 (90..220)	unused							0	FVC	FVC	0	0									0	0			0	FVC	FVC	0	0										0						
Asia	Asia 1		A1	6..90	115..180 (90..220)	unused		0					VC	VC	0	0	0						0	0	0	0	0																							
	Asia 2		A2	6..90	115..180 (90..220)	unused		0					VC	VC	0	0	0						0	0	0	0	0																							
	Asia 3		A3	6..90	115..180 (90..220)	unused		0					VC	VC	0	0	0							0	0	0	0																							
	Asia 4		A4	6..90	115..180 (90..220)	unused		0					VC	VC	0	0	0							0	0	0	0																							
	JRS 2001	2001	JR	18..92 (18..)	90..220	unused		0					0	FVC	0	0	0							0	0	0	0																							

	Reference	Index	Publication
North America	NHANES III (Hankinson)	/REF_P01/	John L. Hankinson, John R. Odencrantz, and Kathleen B. Fedan. Spirometric Reference Values from a Sample of the General U.S. Population. Am J Respir Crit Care Med, Vol 159, p 179-187, 1999
	Knudson_83	/REF_P02/	Knudson, Ronald J, Michael Lebowitz, Holberg Catherine J., Benjamin Burrows. Changes in the Normal Maximal Expiratory Flow-Volume Curve with Aging. American Review of Respiratory Disease, Volume 127, p. 725-734, 1983.
	Knudson_76	/REF_P03/	Knudson, Ronald J, Ronald Slatin, Michael Lebowitz, Benjamin Burrows. The maximal Expiratory Flow-Volume Curve. American Review of Respiratory Disease, Volume 113, p. 587-600, 1976
	Crapo	/REF_P04/	Crapo RO, Morris AH, Gardner RM. Reference spirometric values using techniques and equipment that meets ATS recommendations. Am Rev Respir Dis Volume 123, p.659-664, 1981.
	Morris	/REF_P05/ /REF_P19/	Morris, James F., Koski, Arthur, Lavon Johnson. Spirometric Standards for Healthy Non-Smoking Adults. American Review of Respiratory Disease, Volume 10-3, p. 57-67, 1971 Morris, J.F. West J. Med (1976) 125:110-118
	Hsu	/REF_P06/	Hsu KHK, Bartholomew PH, Thompson V, Hsieh GSJ. Ventilatory Functions of Normal children and Young Adults- Mexican- American, White, Black. I. Spirometry. J Pediatr Volume 95, p. 14-23, 1979.
	Dockery (Harvard)	/REF_P07/	Dockery, D.W. et al., Distribution of Forced Vital Capacity and Forced Expiratory Volume in One Second in Children 6 to 11 Years of Age, American Rev. of Respiratory Disease. Volume 128, p. 405-412, 1983.
	Polgar	/REF_P20/	Polgar, Promadhat, Pulmonary Function Testing in Children: Techniques and Standards. W.B. Saunders Co., Philadelphia, 1971
	Cherniak	/REF_P08/	Cherniak, R.M., and Raber M.B. Normal Standards for Ventilatory Function using an Automated Wedge Spirometer. American Review of Respiratory Disease. Volume 106, p.38-46, 1972
South America	Pereira	/REF_P09/	Carlos Alberto de Castro Pereira, Sueli da Penha Barreto, João Geraldo Simões, Francisco W.L. Pereira, José Gerson Gerstler, Joge Nakatani. Valores de referência para a espirometria em uma amostra da população brasileira adulta, Jornal de Pneumologia 18(1):10-22, maio de 1992
Europe	ERS (ECCS, EGKS)	/REF_P10/	P.H. Quanjer. Lung Volumes and Forced Ventilatory Flows. Eur Respir J, Vol 6, Suppl 16, p. 5-40, 1993
	Zapletal	/REF_P11/	A. Zapletal, T. Paul, M. Samanek. Die Bedeutung heutiger Methoden der Lungenfunktionsdiagnostik zur Feststellung einer Obstruktion der Atemwege bei Kindern und Jugendlichen. Z. Erkrank. Atm.-Org., Volume 149, 343-371, 1977.
	Austria (Forche)	/REF_P12/	G. Forche, K. Harnoncourt, E. Stadlober. Neue spirometrische Bezugswerte für Kinder, Jugendliche und Erwachsene. Öst. Ärztsztg. 43, 15-16, 1988.
	Sapaldia	/REF_P13/	SAPALDIA team, O Brändli, CH. Schindler, N. Künzli, R. Keller, A.P. Perruchoud. Lung function in healthy never smoking adults: reference values and lower limits of normal of a Swiss population. Thorax 1996; 51:277-283
	Spain (Roca)	/REF_P14/	J. Roca et al. spirometric reference values for a Mediterranean population. Bull Eur Physiopathol Respir, 18:101-102, 1982.
Scandinavia	Hedenström	/REF_P15/	H. Hedenström, P. Malmberg, K. Agarwal. Reference values for lung Function tests in females. Bull. Eur. Physiopathol. Respir. 21, p. 551-557, 1985. H. Hedenström, P. Malmberg, H.V. Fridriksson. Reference values for lung function tests in men. Upsala J. Med. Sci., 91:299-310, 1986
	Gulsvik	/REF_P16/	A. Gulsvik. Spirometri (Korrespondanser). Tidsskr Nor Loegeforen nr. 31, 105:2240-2, 1985.
	Berglund	/REF_P17/	E. Berglund, G. Birath, J. Bjure, G. Grimby, I. Kjellmer, L. Sandqvist, B. Söderholm. Spirometric Studies in Normal Subjects. Acta Medica Scandinavica, Vol. 173, fasc. 2, p. 185-206, 1963.
Australasia	Hibbert	/REF_P18/	Marienne E. Hibbert, M App SCI, Anna Lannigan, RN, Louis I. Landau, MD, Peter D. Phelan, MD. Lung Function Values from a Longitudinal Study of Healthy Children Adolescents, Paediatric Pulmonology 7:101-109 (1989)
	Gore, Crockett	/REF_P19/	C.J. Gore, A.J. Crockett, D.G. Pederson, M.L. Booth, A. Bauman, N. Owen. Spirometric standards for healthy adult lifetime nonsmokers in Australia. Eur Respir J., 1995, 8, 773-782
	JRS2001	/REF_P20/	日本人のスパイログラムと動脈血液ガス分圧基準値 日本呼吸器学会肺生理専門委員会 2001年4月

## 20 Appendix D: Electromagnetic Compatibility (EMC)

Changes or modification to the EasyOne system not expressly approved by ndd could cause EMC issues with this or other equipment. The EasyOne system is designed and tested to comply with applicable regulation regarding EMC and needs to be installed and put into service according to the EMC information stated as follows.

### WARNING

Use of portable phones or other radio frequency (RF) emitting equipment near the system may cause unexpected or adverse operation.

### WARNING

The equipment or system should not be used adjacent to, or stacked with, other equipment. If adjacent or stacked use is necessary, the equipment or system should be tested to verify normal operation in the configuration in which it is being used.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The EasyOne is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the EasyOne is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF emissions EN 55011	Group 1	The equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.  The equipment is suitable for use in all establishments inclusively in domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions EN 55011	Class B	
Harmonic Emissions EN 61000-3-2	Not applicable	
Voltage fluctuations/Flicker emissions EN 61000-3-3	Not applicable	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The EasyOne is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to ensure that the EasyOne is used in such an environment.			
Immunity Test	EN 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) EN 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst EN 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	The product has no power supply lines. The product has no input or output lines that require testing.	
Surge EN 61000-4-5	± 1 kV differential mode ± 2 kV common mode	The product has no power supply lines.	
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	<5% $U_t$ (>95% dip in $U_t$ ) for 0.5 cycles 40% $U_t$ (60% dip in $U_t$ ) for 5 cycles 70% $U_t$ (30% dip in $U_t$ ) for 25 cycles <5% $U_t$ (>95% dip in $U_t$ ) for 5 s	The product has no power supply lines	
Power frequency (50/60 Hz) magnetic field EN 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristics of a typical location in a typical commercial or hospital environment.
NOTE: $U_t$ is the AC mains voltage prior to application of the test level.			

**Guidance and Manufacturer's Declaration – Electromagnetic Immunity**

The EasyOne is intended for use in the electromagnetic environment specified below. It is the responsibility of the customer or user to assure that the EasyOne is used in such an environment.

Immunity Test	EN 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Conducted RF EN 61000-4-6 Radiated RF EN 61000-4-3	3 V rms 150 KHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 V rms  3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where <b>P</b> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <b>d</b> is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup>, should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or re-locating the equipment.			
b Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3 V/m.			

**Recommended Separation Distances**

The table below provides the recommended separation distances (in meters) between portable and mobile RF communication equipment and the EasyOne. The EasyOne is intended for use in the electromagnetic environment on which radiated RF disturbances are controlled. The customer or the user of the EasyOne can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the EasyOne as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum Output Power of Transmitter in Watts	Separation Distance in Meters (m) According to Frequency of Transmitter		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance [d] in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note: These guidelines may not apply in all instances. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

**Compliant Cables and Accessories**

The product has no accessories which affect EMC compliance.

ndd Medizintechnik AG  
Technoparkstrasse 1, CH-8005 Zürich, Switzerland  
Tel +41 (44) 445 25 30, Fax +41 (44) 445 25 31, E-Mail info@ndd.ch

## Declaration of Conformity

Manufacturer: ndd Medizintechnik AG  
Address: Technoparkstrasse 1  
CH-8005 Zurich, Switzerland

declares under its sole responsibility, that the product

Product designation: Spirometer

Product name: EasyOne-CS

EasyOne has been classified as Class IIa and is in conformity with the essential requirements and provisions of Council Directive 93/42 EEC,

is in conformity with the following standards transposing harmonized standards:

EN 1041:1998  
EN ISO 14971:2000 and A1 : 2003  
EN 50103:1995  
EN 60601-1: 1990 and A1 :1993 and A2 : 1995  
EN 60601-1-1: 2001  
EN 60601-1-2: 2001  
EN 60601-1-4: 1996 and A1: 1999

is subject to the procedure set out in ISO 9001, ISO 13485, EN 46001 and in Annex 2 of Directive 93/42/EEC under the supervision of Notified Body Number 0120, SGS Yarsley International Certification Services Ltd., 217-221 London Road, Camberley, Surrey, GU15 3EY.

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Copyright© 2007 by nnd Medizintechnik AG, Switzerland. All rights reserved

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	Référence	Index	Publication
Amérique du Nord	NHANES III (Hankinson)	/REF_P01/	John L. Hankinson, John R. Odencrantz, and Kathleen B. Fedan. Spirometric Reference Values from a Sample of the General U.S. Population. Am J Respir Crit Care Med, Vol 159, p 179-187, 1999
	Knudson_83	/REF_P02/	Knudson, Ronald J, Michael Lebowitz, Holberg Catherine J., Benjamin Burrows. Changes in the Normal Maximal Expiratory Flow-Volume Curve with Aging. American Review of Respiratory Disease, Volume 127, p. 725-734, 1983.
	Knudson_76	/REF_P03/	Knudson, Ronald J, Ronald Slatin, Michael Lebowitz, Benjamin Burrows. The maximal Expiratory Flow-Volume Curve. American Review of Respiratory Disease, Volume 113, p. 587-600, 1976
	Crapo	/REF_P04/	Crapo RO, Morris AH, Gardner RM. Reference spirometric values using techniques and equipment that meets ATS recommendations. Am Rev Respir Dis Volume 123, p.659-664, 1981.
	Morris	/REF_P05/ /REF_P19/	Morris, James F., Koski, Arthur, Lavon Johnson. Spirometric Standards for Healthy Non-Smoking Adults. American Review of Respiratory Disease, Volume 10-3, p. 57-67, 1971 Morris, J.F. West J. Med (1976) 125:110-118
	Hsu	/REF_P06/	Hsu KHK, Bartholomew PH, Thompson V, Hsieh GSJ. Ventilatory Functions of Normal children and Young Adults- Mexican- American, White, Black. I. Spirometry. J Pediatr Volume 95, p. 14-23, 1979.
	Dockery (Harvard)	/REF_P07/	Dockery, D.W. et al., Distribution of Forced Vital Capacity and Forced Expiratory Volume in One Second in Children 6 to 11 Years of Age, American Rev. of Respiratory Disease. Volume 128, p. 405-412, 1983.
	Polgar	/REF_P20/	Polgar, Promadhat, Pulmonary Function Testing in Children: Techniques and Standards. W.B. Saunders Co., Philadelphia, 1971
	Cherniak	/REF_P08/	Cherniak, R.M., and Raber M.B. Normal Standards for Ventilatory Function using an Automated Wedge Spirometer. American Review of Respiratory Disease. Volume 106, p.38-46, 1972
Amérique du Sud	Pereira	/REF_P09/	Carlos Alberto de Castro Pereira, Sueli da Penha Barreto, João Geraldo Simões, Francisco W.L. Pereira, José Gerson Gerstler, Joge Nakatani. Valores de referência para a espirometria em uma amostra da população brasileira adulta, Jornal de Pneumologia 18(1):10-22, maio de 1992
Europe	ERS (ECCS, EGKS)	/REF_P10/	P.H. Quanjer. Lung Volumes and Forced Ventilatory Flows. Eur Respir J, Vol 6, Suppl 16, p. 5-40, 1993
	Zapletal	/REF_P11/	A. Zapletal, T. Paul, M. Samanek. Die Bedeutung heutiger Methoden der Lungenfunktionsdiagnostik zur Feststellung einer Obstruktion der Atemwege bei Kindern und Jugendlichen. Z. Erkrank. Atm.-Org., Volume 149, 343-371, 1977.
	Austria (Forche)	/REF_P12/	G. Forche, K. Harnoncourt, E. Stadlober. Neue spirometrische Bezugswerte für Kinder, Jugendliche und Erwachsene. Öst. Ärztztg. 43, 15-16, 1988.
	Sapaldia	/REF_P13/	SAPALDIA team, O Brändli, CH. Schindler, N. Künzli, R. Keller, A.P. Perruchoud. Lung function in healthy never smoking adults: reference values and lower limits of normal of a Swiss population. Thorax 1996; 51:277-283
	Spain (Roca)	/REF_P14/	J. Roca et al. spirometric reference values for a Mediterranean population. Bull Eur Physiopathol Respir, 18:101-102, 1982.
Scandinavie	Hedenström	/REF_P15/	H. Hedenström, P. Malmberg, K. Agarwal. Reference values for lung Function tests in females. Bull. Eur. Physiopathol. Respir. 21, p. 551-557, 1985. H. Hedenström, P. Malmberg, H.V. Fridriksson. Reference values for lung function tests in men. Upsala J. Med. Sci., 91:299-310, 1986
	Gulsvik	/REF_P16/	A. Gulsvik. Spirometri (Korrespondanser). Tidsskr Nor Loegeforen nr. 31, 105:2240-2, 1985.
	Berglund	/REF_P17/	E. Berglund, G. Birath, J. Bjure, G. Grimby, I. Kjellmer, L. Sandqvist, B. Söderholm. Spirometric Studies in Normal Subjects. Acta Medica Scandinavica, Vol. 173, fasc. 2, p. 185-206, 1963.
Australie Asie	Hibbert	/REF_P18/	Marienne E. Hibbert, M App SCI, Anna Lannigan, RN, Louis I. Landau, MD, Peter D. Phelan, MD. Lung Function Values From a Longitudinal Study of Healthy Children Adolescents, Pediatric Pulmonology 7:101-109 (1989)
	Gore, Crockett	/REF_P19/	C.J. Gore, A.J. Crockett, D.G. Pederson, M.L. Booth, A. Bauman, N. Owen. Spirometric standards for healthy adult lifetime nonsmokers in Australia. Eur Respir J., 1995, 8, 773-782
	JRS2001	/REF_P20/	日本人のスパイログラムと動脈血液ガス分圧基準値 日本呼吸器学会肺生理専門委員会 2001年4月

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	Riferimento	Indice	Pubblicazione
America del nord	NHANES III (Hankinson)	/REF_P01/	John L. Hankinson, John R. Odencrantz, and Kathleen B. Fedan. Spirometric Reference Values from a Sample of the General U.S. Population. Am J Respir Crit Care Med, Vol 159, p 179-187, 1999
	Knudson_83	/REF_P02/	Knudson, Ronald J, Michael Lebowitz, Holberg Catherine J., Benjamin Burrows. Changes in the Normal Maximal Expiratory Flow-Volume Curve with Aging. American Review of Respiratory Disease, Volume 127, p. 725-734, 1983.
	Knudson_76	/REF_P03/	Knudson, Ronald J, Ronald Slatin, Michael Lebowitz, Benjamin Burrows. The maximal Expiratory Flow-Volume Curve. American Review of Respiratory Disease, Volume 113, p. 587-600, 1976
	Crapo	/REF_P04/	Crapo RO, Morris AH, Gardner RM. Reference spirometric values using techniques and equipment that meets ATS recommendations. Am Rev Respir Dis Volume 123, p.659-664, 1981.
	Morris	/REF_P05/ /REF_P19/	Morris, James F., Koski, Arthur, Lavon Johnson. Spirometric Standards for Healthy Non-Smoking Adults. American Review of Respiratory Disease, Volume 10-3, p. 57-67, 1971 Morris, J.F. West J. Med (1976) 125:110-118
	Hsu	/REF_P06/	Hsu KHK, Bartholomew PH, Thompson V, Hsieh GSJ. Ventilatory Functions of Normal children and Young Adults- Mexican- American, White, Black. I. Spirometry. J Pediatr Volume 95, p. 14-23, 1979.
	Dockery (Harvard)	/REF_P07/	Dockery, D.W. et al., Distribution of Forced Vital Capacity and Forced Expiratory Volume in One Second in Children 6 to 11 Years of Age, American Rev. of Respiratory Disease. Volume 128, p. 405-412, 1983.
	Polgar	/REF_P20/	Polgar, Promadhat, Pulmonary Function Testing in Children: Techniques and Standards. W.B. Saunders Co., Philadelphia, 1971
	Cherniak	/REF_P08/	Cherniak, R.M., and Raber M.B. Normal Standards for Ventilatory Function using an Automated Wedge Spirometer. American Review of Respiratory Disease. Volume 106, p.38-46, 1972
America del sud	Pereira	/REF_P09/	Carlos Alberto de Castro Pereira, Sueli da Penha Barreto, João Geraldo Simões, Francisco W.L. Pereira, José Gerson Gerstler, Joge Nakatani. Valores de referência para a espirometria em uma amostra da população brasileira adulta, Jornal de Pneumologia 18(1):10-22, maio de 1992
Europa	ERS (ECCS, EGKS)	/REF_P10/	P.H. Quanjer. Lung Volumes and Forced Ventilatory Flows. Eur Respir J, Vol 6, Suppl 16, p. 5-40, 1993
	Zapletal	/REF_P11/	A. Zapletal, T. Paul, M. Samanek. Die Bedeutung heutiger Methoden der Lungenfunktionsdiagnostik zur Feststellung einer Obstruktion der Atemwege bei Kindern und Jugendlichen. Z. Erkrank. Atm.-Org., Volume 149, 343-371, 1977.
	Austria (Forche)	/REF_P12/	G. Forche, K. Harnoncourt, E. Stadlober. Neue spirometrische Bezugswerte für Kinder, Jugendliche und Erwachsene. Öst. Ärztsztg. 43, 15-16, 1988.
	Sapaldia	/REF_P13/	SAPALDIA team, O Brändli, CH. Schindler, N. Künzli, R. Keller, A.P. Perruchoud. Lung function in healthy never smoking adults: reference values and lower limits of normal of a Swiss population. Thorax 1996; 51:277-283
	Spain (Roca)	/REF_P14/	J. Roca et al. spirometric reference values for a Mediterranean population. Bull Eur Physiopathol Respir, 18:101-102, 1982.
Scandinavia	Hedenström	/REF_P15/	H. Hedenström, P. Malmberg, K. Agarwal. Reference values for lung Function tests in females. Bull. Eur. Physiopathol. Respir. 21, p. 551-557, 1985. H. Hedenström, P. Malmberg, H.V. Fridriksson. Reference values for lung function tests in men. Upsala J. Med. Sci., 91:299-310, 1986
	Gulsvik	/REF_P16/	A. Gulsvik. Spirometri (Korrespondanser). Tidsskr Nor Loegeforen nr. 31, 105:2240-2, 1985.
	Berglund	/REF_P17/	E. Berglund, G. Birath, J. Bjure, G. Grimby, I. Kjellmer, L. Sandqvist, B. Söderholm. Spirometric Studies in Normal Subjects. Acta Medica Scandinavica, Vol. 173, fasc. 2, p. 185-206, 1963.
Australia/Asia	Hibbert	/REF_P18/	Marionne E. Hibbert, M App SCI, Anna Lannigan, RN, Louis I. Landau, MD, Peter D. Phelan, MD. Lung Function Values From a Longitudinal Study of Healthy Children Adolescents, Pediatric Pulmonology 7:101-109 (1989)
	Gore, Crockett	/REF_P19/	C.J. Gore, A.J. Crockett, D.G. Pederson, M.L. Booth, A. Bauman, N. Owen. Spirometric standards for healthy adult lifetime nonsmokers in Australia. Eur Respir J., 1995, 8, 773-782
	JRS2001	/REF_P20/	日本人のスパイログラムと動脈血液ガス分圧基準値 日本呼吸器学会肺生理専門委員会 2001年4月

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	Referencia	Índice	Publicación
América del Norte	NHANES III (Hankinson)	/REF_P01/	John L. Hankinson, John R. Odencrantz, and Kathleen B. Fedan. Spirometric Reference Values from a Sample of the General U.S. Population. Am J Respir Crit Care Med, Vol 159, p 179-187, 1999
	Knudson_83	/REF_P02/	Knudson, Ronald J, Michael Lebowitz, Holberg Catherine J., Benjamin Burrows. Changes in the Normal Maximal Expiratory Flow-Volume Curve with Aging. American Review of Respiratory Disease, Volume 127, p. 725-734, 1983.
	Knudson_76	/REF_P03/	Knudson, Ronald J, Ronald Slatin, Michael Lebowitz, Benjamin Burrows. The maximal Expiratory Flow-Volume Curve. American Review of Respiratory Disease, Volume 113, p. 587-600, 1976
	Crapo	/REF_P04/	Crapo RO, Morris AH, Gardner RM. Reference spirometric values using techniques and equipment that meets ATS recommendations. Am Rev Respir Dis Volume 123, p.659-664, 1981.
	Morris	/REF_P05/ /REF_P19/	Morris, James F., Koski, Arthur, Lavon Johnson. Spirometric Standards for Healthy Non-Smoking Adults. American Review of Respiratory Disease, Volume 10-3, p. 57-67, 1971 Morris, J.F. West J. Med (1976) 125:110-118
	Hsu	/REF_P06/	Hsu KHK, Bartholomew PH, Thompson V, Hsieh GSJ. Ventilatory Functions of Normal children and Young Adults- Mexican- American, White, Black. I. Spirometry. J Pediatr Volume 95, p. 14-23, 1979.
	Dockery (Harvard)	/REF_P07/	Dockery, D.W. et al., Distribution of Forced Vital Capacity and Forced Expiratory Volume in One Second in Children 6 to 11 Years of Age, American Rev. of Respiratory Disease. Volume 128, p. 405-412, 1983.
	Polgar	/REF_P20/	Polgar, Promadhat, Pulmonary Function Testing in Children: Techniques and Standards. W.B. Saunders Co., Philadelphia, 1971
	Cherniak	/REF_P08/	Cherniak, R.M., and Raber M.B. Normal Standards for Ventilatory Function using an Automated Wedge Spirometer. American Review of Respiratory Disease. Volume 106, p.38-46, 1972
América del Sur	Pereira	/REF_P09/	Carlos Alberto de Castro Pereira, Sueli da Penha Barreto, João Geraldo Simões, Francisco W.L. Pereira, José Gerson Gerstler, Joge Nakatani. Valores de referência para a espirometria em uma amostra da população brasileira adulta, Jornal de Pneumologia 18(1):10-22, maio de 1992
Europa	ERS (ECCS, EGKS)	/REF_P10/	P.H. Quanjer. Lung Volumes and Forced Ventilatory Flows. Eur Respir J, Vol 6, Suppl 16, p. 5-40, 1993
	Zapletal	/REF_P11/	A. Zapletal, T. Paul, M. Samanek. Die Bedeutung heutiger Methoden der Lungenfunktionsdiagnostik zur Feststellung einer Obstruktion der Atemwege bei Kindern und Jugendlichen. Z. Erkrank. Atm.-Org., Volume 149, 343-371, 1977.
	Austria (Forche)	/REF_P12/	G. Forche, K. Harnoncourt, E. Stadlober. Neue spirometrische Bezugswerte für Kinder, Jugendliche und Erwachsene. Öst. Ärztztg. 43, 15-16, 1988.
	Sapaldia	/REF_P13/	SAPALDIA team, O Brändli, CH. Schindler, N. Künzli, R. Keller, A.P. Perruchoud. Lung function in healthy never smoking adults: reference values and lower limits of normal of a Swiss population. Thorax 1996; 51:277-283
	Spain (Roca)	/REF_P14/	J. Roca et al. spirometric reference values for a Mediterranean population. Bull Eur Physiopathol Respir, 18:101-102, 1982.
Escandinavia	Hedenström	/REF_P15/	H. Hedenström, P. Malmberg, K. Agarwal. Reference values for lung Function tests in females. Bull. Eur. Physiopathol. Respir. 21, p. 551-557, 1985. H. Hedenström, P. Malmberg, H.V. Fridriksson. Reference values for lung function tests in men. Upsala J. Med. Sci., 91:299-310, 1986
	Gulsvik	/REF_P16/	A. Gulsvik. Spirometri (Korrespondanser). Tidsskr Nor Loegeforen nr. 31, 105:2240-2, 1985.
	Berglund	/REF_P17/	E. Berglund, G. Birath, J. Bjure, G. Grimby, I. Kjellmer, L. Sandqvist, B. Söderholm. Spirometric Studies in Normal Subjects. Acta Medica Scandinavica, Vol. 173, fasc. 2, p. 185-206, 1963.
Australia - Asia	Hibbert	/REF_P18/	Marianne E. Hibbert, M App SCI, Anna Lannigan, RN, Louis I. Landau, MD, Peter D. Phelan, MD. Lung Function Values From a Longitudinal Study of Healthy Children Adolescents, Pediatric Pulmonology 7:101-109 (1989)
	Gore, Crockett	/REF_P19/	C.J. Gore, A.J. Crockett, D.G. Pederson, M.L. Booth, A. Bauman, N. Owen. Spirometric standards for healthy adult lifetime nonsmokers in Australia. Eur Respir J., 1995, 8, 773-782
	JRS2001	/REF_P20/	日本人のスパイログラムと動脈血液ガス分圧基準値 日本呼吸器学会肺生理専門委員会 2001年4月

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	Reference	Index	Publication
North America	NHANES III (Hankinson)	/REF_P01/	John L. Hankinson, John R. Odencrantz, and Kathleen B. Fedan. Spirometric Reference Values from a Sample of the General U.S. Population. Am J Respir Crit Care Med, Vol 159, p 179-187, 1999
	Knudson_83	/REF_P02/	Knudson, Ronald J, Michael Lebowitz, Holberg Catherine J., Benjamin Burrows. Changes in the Normal Maximal Expiratory Flow-Volume Curve with Aging. American Review of Respiratory Disease, Volume 127, p. 725-734, 1983.
	Knudson_76	/REF_P03/	Knudson, Ronald J, Ronald Slatin, Michael Lebowitz, Benjamin Burrows. The maximal Expiratory Flow-Volume Curve. American Review of Respiratory Disease, Volume 113, p. 587-600, 1976
	Crapo	/REF_P04/	Crapo RO, Morris AH, Gardner RM. Reference spirometric values using techniques and equipment that meets ATS recommendations. Am Rev Respir Dis Volume 123, p.659-664, 1981.
	Morris	/REF_P05/ /REF_P19/	Morris, James F., Koski, Arthur, Lavon Johnson. Spirometric Standards for Healthy Non-Smoking Adults. American Review of Respiratory Disease, Volume 10-3, p. 57-67, 1971 Morris, J.F. West J. Med (1976) 125:110-118
	Hsu	/REF_P06/	Hsu KHK, Bartholomew PH, Thompson V, Hsieh GSJ. Ventilatory Functions of Normal children and Young Adults- Mexican- American, White, Black. I. Spirometry. J Pediatr Volume 95, p. 14-23, 1979.
	Dockery (Harvard)	/REF_P07/	Dockery, D.W. et al., Distribution of Forced Vital Capacity and Forced Expiratory Volume in One Second in Children 6 to 11 Years of Age, American Rev. of Respiratory Disease. Volume 128, p. 405-412, 1983.
	Polgar	/REF_P20/	Polgar, Promadhat, Pulmonary Function Testing in Children: Techniques and Standards. W.B. Saunders Co., Philadelphia, 1971
	Cherniak	/REF_P08/	Cherniak, R.M., and Raber M.B. Normal Standards for Ventilatory Function using an Automated Wedge Spirometer. American Review of Respiratory Disease. Volume 106, p.38-46, 1972
South America	Pereira	/REF_P09/	Carlos Alberto de Castro Pereira, Sueli da Penha Barreto, João Geraldo Simões, Francisco W.L. Pereira, José Gerson Gerstler, Joge Nakatani. Valores de referência para a espirometria em uma amostra da população brasileira adulta, Jornal de Pneumologia 18(1):10-22, maio de 1992
Europe	ERS (ECCS, EGKS)	/REF_P10/	P.H. Quanjer. Lung Volumes and Forced Ventilatory Flows. Eur Respir J, Vol 6, Suppl 16, p. 5-40, 1993
	Zapletal	/REF_P11/	A. Zapletal, T. Paul, M. Samanek. Die Bedeutung heutiger Methoden der Lungenfunktionsdiagnostik zur Feststellung einer Obstruktion der Atemwege bei Kindern und Jugendlichen. Z. Erkrank. Atm.-Org., Volume 149, 343-371, 1977.
	Austria (Forche)	/REF_P12/	G. Forche, K. Harnoncourt, E. Stadlober. Neue spirometrische Bezugswerte für Kinder, Jugendliche und Erwachsene. Öst. Ärztsztg. 43, 15-16, 1988.
	Sapaldia	/REF_P13/	SAPALDIA team, O Brändli, CH. Schindler, N. Künzli, R. Keller, A.P. Perruchoud. Lung function in healthy never smoking adults: reference values and lower limits of normal of a Swiss population. Thorax 1996; 51:277-283
	Spain (Roca)	/REF_P14/	J. Roca et al. spirometric reference values for a Mediterranean population. Bull Eur Physiopathol Respir, 18:101-102, 1982.
Scandinavia	Hedenström	/REF_P15/	H. Hedenström, P. Malmberg, K. Agarwal. Reference values for lung Function tests in females. Bull. Eur. Physiopathol. Respir. 21, p. 551-557, 1985. H. Hedenström, P. Malmberg, H.V. Fridriksson. Reference values for lung function tests in men. Upsala J. Med. Sci., 91:299-310, 1986
	Gulsvik	/REF_P16/	A. Gulsvik. Spirometri (Korrespondanser). Tidsskr Nor Loegeforen nr. 31, 105:2240-2, 1985.
	Berglund	/REF_P17/	E. Berglund, G. Birath, J. Bjure, G. Grimby, I. Kjellmer, L. Sandqvist, B. Söderholm. Spirometric Studies in Normal Subjects. Acta Medica Scandinavica, Vol. 173, fasc. 2, p. 185-206, 1963.
Australasia	Hibbert	/REF_P18/	Marionne E. Hibbert, M App SCI, Anna Lannigan, RN, Louis I. Landau, MD, Peter D. Phelan, MD. Lung Function Values From a Longitudinal Study of Healthy Children Adolescents, Paediatric Pulmonology 7:101-109 (1989)
	Gore, Crockett	/REF_P19/	C.J. Gore, A.J. Crockett, D.G. Pederson, M.L. Booth, A. Bauman, N. Owen. Spirometric standards for healthy adult lifetime nonsmokers in Australia. Eur Respir J., 1995, 8, 773-782
	JRS2001	/REF_P20/	日本人のスパイログラムと動脈血液ガス分圧基準値 日本呼吸器学会肺生理専門委員会 2001年4月

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---

*European Headquarters  
GE Medical Systems  
Information Technologies GmbH  
Munzinger Straße 3  
79111 Freiburg  
Germany  
Tel. +49 761 45 43 – 0  
Fax +49 761 45 43 – 233*

*World Headquarters  
GE Medical Systems  
Information Technologies,  
8200 West Tower Avenue  
Milwaukee, WI 53223  
United States of America  
Tel.: + 1 414 355 5000  
Fax +1 414 355 3790*

*Asian Headquarters  
GE Medical Systems  
Inc. Information Technologies Asia  
GE (China) Co., Ltd.  
24th Floor Shanghai MAXDO Centre,  
8 Xing Yi Road, Hong Qiao Development Zone  
Shanghai 200336, People's Republic of China  
Tel: +86-21-5257-4650  
Fax: +86-21-5208-2008*

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Copyright© 2007 by ndd Medizintechnik AG, Switzerland. Tutti i diritti riservati
Copyright© 2007 by nnd Medizintechnik AG, Suiza. Reservados todos los derechos.
Copyright© 2007 by ndd Medizintechnik AG, Suíça. Todos os direitos reservados